

Multi-Lumen CVC Kit with Blue FlexTip® Catheter

Product No.

ASK-12703-PRJ



Contents:

- 1: Multi-Lumen Indwelling Catheter: 7 Fr. x 16 cm Radiopaque Polyurethane with Blue FlexTip®, Extension Line Clamps
- 1: Spring-Wire Guide, Marked: .032" (0.81 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer™
- 1: Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle
- 1: Pressure Transduction Probe
- 1: Syringe: 10 mL Luer-Lock
- 1: Introducer Needle: Echogenic 18 Ga. x 2-1/2" (6.35 cm) XTW and 5 mL Arrow® Raulerson Spring-Wire Introduction Syringe
- 1: Injection Needle: Eclipse™¹ 22 Ga. x 1-1/2" (3.81 cm) RW and 5 mL Luer-Slip Syringe
- 1: Injection Needle: Eclipse™¹ 25 Ga. x 1" (2.54 cm) RW and 3 mL Luer-Lock Syringe
- 1: Tissue Dilator
- 1: Needle Holder
- 1: Checklist/CLIP Sheet
- 1: Sterile Procedure Sign
- 1: Maximal Barrier Drape™ with 4" fenestration
- 1: CSR Wrap
- 1: Forceps
- 1: Suture: 3-0 Silk with Curved Needle
- 2: Gauze Pads: 2" x 2"
- 5: Gauze Pads: 4" x 4"
- 1: SharpsAway® Disposal Cup
- 1: SecondSite™ Adjustable Hub: Catheter Clamp and Fastener
- 1: Safety Scalpel: #11
- 1: SharpsAway II™ Locking Disposal Cup
- 1: Dressing: Tegaderm®² 10 cm x 12 cm
- 1: Towel
- 1: HemoHopper® Fluid Receptacle


¹ A trademark of Becton, Dickinson and Company.
² A registered trademark of 3M Company.

Arrow Select Kits



Rx only. Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.

Sterile. Single use – do not resterilize. Sterilized by ethylene oxide. Keep away from sunlight. Keep dry. Do not use if package has been previously opened or damaged. Not made with natural rubber latex. Fluid path components are non-pyrogenic. Contains no medication.




Lumen	Priming Volume* (mL)	Flow Rate† (mL/hr)
Distal (16 Ga.)	0.40	3600
Medial (18 Ga.)	0.36	1850
Proximal (18 Ga.)	0.40	2020

* Priming volumes are approximate and are done without an injection cap.
 † Based on 100 cm head height.

Arrow International, Inc. | 2400 Bernville Road, Reading, PA 19605 USA | www.arrowintl.com

J-12703-102B (9/16)



Multi-Lumen CVC Kit with Blue FlexTip® Catheter

Product No.

ASK-12703-PRJ



7 Fr.

3 Lumen

16 cm catheter length

.032 inch dia. spring-wire guide

