



REF CV-12122-FX

(IPNXXXXXX)

12 Fr. | 2 Lumen | 16 cm catheter length | .035 inch dia. spring-wire guide

Two-Lumen Hemodialysis Catheterization Set with Blue FlexTip® Catheter for High Volume Infusions

Contents:


- 1: Two-Lumen Indwelling Catheter: 12 Fr. x 16 cm Radiopaque Polyurethane with Blue FlexTip®, Extension Line Clamps
- 1: Spring-Wire Guide, Marked: .035" (0.89 mm) dia. x 23-5/8" (60 cm) (Straight Soft Tip on One End – "J" Tip on Other) with Arrow Advancer
- 1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW
- 1: Syringe: 5 mL Luer-Slip
- 1: Tissue Dilator: 12 Fr. (4.0 mm) x 14 cm
- 2: Dust Cap: Non-Vented

Rx only. Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.

Not made with natural rubber latex.

Fluid path components are non-pyrogenic.

Contraindications: The Arrow® large-bore two-lumen catheter is not designed for long-term hemodialysis or for use in patients with thrombosed vessels.



| Lumen | Priming Volume* (ml) | Flow Rate† (ml/min) | Venous Pressure Setting (mmHg) |
|-------------------|----------------------|---------------------|--------------------------------|
| Distal (12 Ga.) | 1.4 | 220 | 100 |
| | | 370 | 200 |
| Proximal (12 Ga.) | 1.3 | 220 | 100 |
| | | 370 | 200 |

* Priming volumes are approximate and are done without accessories.

† Blood flow rates are approximate and were determined by varying the venous pressure setting from 100 to 200 mmHg in an in vitro experiment using blood analog solution.



LOT LotNumber
YYYY-MM-DD



(01)XXXXXXXXXXXXXX
(17)YYMMDD
(10)LotNumber

Arrow International, Inc.
Subsidiary of Teleflex Incorporated
2400 Bernville Road
Reading, PA 19605 USA

Teleflex®
www.teleflex.com

LBL010059 (11/18)



REF CV-12122-FX 12 Fr. | 2 Lumen | 16 cm catheter length | .035 inch dia. spring-wire guide



Two-Lumen Hemodialysis Catheterization Set with Blue FlexTip® Catheter for High Volume Infusions