



REF BR-15854 (IPNXXXXXX)

8.5^{Fr.} **4^{Lumen}** **20^{cm}** catheter length **.032^{inch dia.}** spring-wire guide

Quad-Lumen Central Venous Catheterization Set with Blue FlexTip® Catheter

Contents:

- 1: Quad-Lumen Indwelling Catheter: 8.5 Fr. x 20 cm Radiopaque Polyurethane with Blue FlexTip®, Extension Line Clamps
- 1: Spring-Wire Guide, Marked: .032" (0.81 mm) dia. x 23-5/8" (60 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer with ECG Mark
- 1: Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle
- 1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW
- 1: Syringe: 5 ml Luer-Slip
- 1: Tissue Dilator: 9 Fr. (3.0 mm) x 10.2 cm
- 4: Dust Cap: Non-Vented
- 1: SecondSite™ Adjustable Hub: Fastener
- 1: SecondSite™ Adjustable Hub: Catheter Clamp
- 1: Scalpel: #11
- 4: Stopcock: 4-way High-Flow


- fr** Jeu de cathétérisme veineux central à quatre lumières avec cathéter Blue FlexTip
- de** Besteck für vierlumigen, zentralen Venenkatheter mit Blue-FlexTip-Katheter
- it** Set per cateterismo della vena centrale a quattro lumi con catetere Blue FlexTip
- pt** Conjunto de Cateterização Venosa Central com Lúmen quádruplo e Cateter Blue FlexTip
- es** Conjunto de cateterización venosa central de cuatro luces con catéter Blue FlexTip
- sv** Set för fyrkanalig central venkatetrering med Blue FlexTip-kateter

All components are CE 0086 unless otherwise noted.

Rx only. Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.

Not made with natural rubber latex.

Fluid path components are non-pyrogenic.



Lumen	Priming Volume* (ml)	Flow Rate† (ml/hr)
Distal (16 Ga.)	0.43	2700
Medial 1 (14 Ga.)	0.54	5100
Medial 2 (18 Ga.)	0.36	1400
Proximal (18 Ga.)	0.41	1600

* Priming volumes are approximate and are done without accessories.

† Based on 100 cm head height.



LOT LotNumber
YYYY-MM-DD



(01)XXXXXXXXXXXXXX
(17)YYMMDD
(10)LotNumber

EC REP Teleflex Medical
IDA Business and
Technology Park
Dublin Road, Athlone
Co. Westmeath, Ireland

Arrow International, Inc.
Subsidiary of Teleflex Incorporated
2400 Bernville Road
Reading, PA 19605 USA

Teleflex®
www.teleflex.com

LBL010372 (01/19)



REF BR-15854 **8.5^{Fr.}** **4^{Lumen}** **20^{cm}** catheter length **.032^{inch dia.}** spring-wire guide

Quad-Lumen Central Venous Catheterization Set with Blue FlexTip® Catheter

