



**REF CS-10853** (IPNXXXXXX)

**8.5** Fr. **3** Lumen **11** cm catheter length **.032** inch dia. spring-wire guide

# Multi-Lumen Central Venous Catheterization Set with Blue FlexTip® Catheter

**Contents:**

- 1: Multi-Lumen Indwelling Catheter: 8.5 Fr. x 11 cm Radiopaque Polyurethane with Blue FlexTip®, Extension Line Clamps
- 1: Spring-Wire Guide, Marked: .032" (0.81 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End – "J" Tip on Other) with Arrow Advancer
- 1: Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle
- 1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW
- 1: Syringe: 5 mL Luer-Slip
- 1: Tissue Dilator: 9 Fr. (3.0 mm) x 10.2 cm
- 3: Dust Cap: Non-Vented
- 1: SecondSite™ Adjustable Hub: Fastener
- 1: SecondSite™ Adjustable Hub: Catheter Clamp

All components are CE 0086 unless otherwise noted.

**Rx only. Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.**


**California Prop. 65**

**⚠ WARNING:** Cancer and Reproductive Harm. [www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov)

Not made with natural rubber latex.

Fluid path components are non-pyrogenic.

- cs** Vicelumenový centrální žilní katetizační set s katetrem Blue FlexTip
- da** Multilumen centralvenekateriseringssett med Blue FlexTip kateter
- nl** Multilumen centraalveneuze katheterisatieset met Blue FlexTip katheter
- fr** Jeu de cathétérisme veineux central à lumières multiples avec cathéter Blue FlexTip
- de** Besteck für viellumigen zentralen Venenkatheter mit Blue FlexTip Katheter
- el** Σετ κεντρικού φλεβικού καθετήριασμού πολλαπλών αυλών με καθετήρα Blue FlexTip
- it** Set per cateterismo della vena centrale multilume con catetere Blue FlexTip
- no** Multilumen sentralt venekateriseringssett med Blue FlexTip kateter
- pl** Zestaw wielokanalowy do cewnikowania żył centralnych z cewnikiem typu Blue FlexTip
- pt** Conjunto de cateterização venosa central com lúmen múltiplo e cateter Blue FlexTip
- ru** Набор с многопросветным центральным венозным катетером с кончиком Blue FlexTip
- sk** Multilúmenová centrálna venózna katetizačná sada s katétrom Blue FlexTip
- sl** Večlumenski centralni venski komplet za katetizacijo s katetrom Blue FlexTip
- es** Conjunto de cateterización venosa central de luz múltiple con catéter Blue FlexTip
- sv** Set för flerkanalig central venkatetrisering med Blue FlexTip kateter
- tr** Multi-Lümen Santral Venöz Kateterizasyon Seti, Blue FlexTip Kateter ile



| Lumen             | Priming Volume* (ml) | Flow Rate† (ml/hr) |
|-------------------|----------------------|--------------------|
| Distal (16 Ga.)   | 0.37                 | 4235               |
| Medial (14 Ga.)   | 0.72                 | 9700               |
| Proximal (16 Ga.) | 0.39                 | 4300               |

\* Priming volumes are approximate and are done without accessories.

† Based on 100 cm head height.



**LOT** LotNumber  
YYYY-MM-DD



(01)XXXXXXXXXXXXXX  
(17)YYMMDD  
(10)LotNumber

**EC REP** Teleflex Medical  
IDA Business and  
Technology Park  
Dublin Road, Athlone  
Co. Westmeath, Ireland

Arrow International, Inc.  
Subsidiary of Teleflex Incorporated  
2400 Bernville Road  
Reading, PA 19605 USA

**Teleflex®**  
[www.teleflex.com](http://www.teleflex.com)

LBL010377 (01/19)



**REF CS-10853** **8.5** Fr. **3** Lumen **11** cm catheter length **.032** inch dia. spring-wire guide

**Multi-Lumen Central Venous Catheterization Set with Blue FlexTip® Catheter**

