

**ARROW**



**REF CS-12854-E** (IPNXXXXXX)

**8.5** Fr. | **4** Lumen | **16** cm catheter length | **.032** inch dia. spring-wire guide

# Quad-Lumen Central Venous Catheterization Set with Blue FlexTip® Catheter

**Contents:**

- 1: Quad-Lumen Indwelling Catheter: 8.5 Fr. x 16 cm Radiopaque Polyurethane with Blue FlexTip®, Extension Line Clamps
- 1: Spring-Wire Guide, Marked: .032" (0.81 mm) dia. x 23-5/8"(60 cm) (Straight Soft Tip on One End – "J" Tip on Other) with Arrow Advancer with ECG Mark
- 1: Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle
- 1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW
- 1: Pressure Transduction Probe
- 1: Syringe: 5 mL Luer-Slip
- 1: Arrow® Raulerson Spring-Wire Introduction Syringe: 5 mL
- 1: Tissue Dilator: 9 Fr. (3.0 mm) x 10.2 cm
- 4: Dust Cap: Non-Vented
- 1: SecondSite™ Adjustable Hub: Fastener
- 1: SecondSite™ Adjustable Hub: Catheter Clamp


- fr** Jeu de cathétérisme veineux central à quatre lumières avec cathéter Blue FlexTip
- de** Besteck für vierlumigen, zentralen Venenkatheter mit Blue FlexTip Katheter
- it** Set per cateterismo della vena centrale a quattro lumi con catetere Blue FlexTip
- pt** Conjunto de Cateterização Venosa Central com Lúmen quádruplo e Cateter Blue FlexTip
- es** Conjunto de cateterización venosa central de cuatro luces con catéter Blue FlexTip
- sv** Set för fyrkanalig central venkatetrering med Blue FlexTip kateter

All components are CE 0086 unless otherwise noted.

**Rx only. Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.**

Not made with natural rubber latex.

Fluid path components are non-pyrogenic.



Lumen	Priming Volume* (ml)	Flow Rate† (ml/hr)
Distal (16 Ga.)	0.39	3200
Medial 1 (14 Ga.)	0.47	5400
Medial 2 (18 Ga.)	0.36	1600
Proximal (18 Ga.)	0.38	1800

\* Priming volumes are approximate and are done without accessories.

† Based on 100 cm head height.



**LOT** LotNumber  
YYYY-MM-DD



(01)XXXXXXXXXXXXXX  
(17)YYMMDD  
(10)LotNumber

**EC REP** Teleflex Medical  
IDA Business and  
Technology Park  
Dublin Road, Athlone  
Co. Westmeath, Ireland

Arrow International, Inc.  
Subsidiary of Teleflex Incorporated  
2400 Bernville Road  
Reading, PA 19605 USA

**Teleflex**  
www.teleflex.com

LBL010385 (11/18)



**REF CS-12854-E** **8.5** Fr. | **4** Lumen | **16** cm catheter length | **.032** inch dia. spring-wire guide

**Quad-Lumen Central Venous Catheterization Set with Blue FlexTip® Catheter**

**ARROW**