



REF CS-14854-E (IPNXXXXXX)

8.5 Fr. 4 Lumen 30 cm catheter length .032 inch dia. spring-wire guide

Quad-Lumen Central Venous Catheterization Set with Blue FlexTip® Catheter

Contents:

- 1: Four-Lumen Catheter: 8.5 Fr. (3.0 mm OD) x 30 cm
- 1: Spring-Wire Guide, Marked: .032" (0.81 mm) dia. x 26-3/4" (68 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer with ECG Mark
- 1: Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle
- 1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW
- 1: Pressure Transduction Probe
- 1: Syringe: 5 mL Luer-Slip
- 1: Arrow® Raulerson Spring-Wire Introduction Syringe: 5 mL
- 1: Tissue Dilator: 10 Fr. (3.5 mm) x 10.2 cm
- 4: Dust Cap: Non-Vented
- 1: SecondSite™ Adjustable Hub: Fastener
- 1: SecondSite™ Adjustable Hub: Catheter Clamp

All components are CE 0086 unless otherwise noted.

- fr** Jeu de cathétérisme veineux central à quatre lumières avec cathéter Blue FlexTip
- de** Besteck für vierlumigen, zentralen Venenkatheter mit Blue FlexTip Katheter
- it** Set per cateterismo della vena centrale a quattro lumi con catetere Blue FlexTip
- pt** Conjunto de Cateterização Venosa Central com Lúmen Quádruplo e Cateter Blue FlexTip
- es** Conjunto de cateterización venosa central de cuatro luces con catéter Blue FlexTip
- sv** Set för fyrkanalig central venkatetrisering med Blue FlexTip-kateter

Rx only. Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.

Not made with natural rubber latex.

Fluid path components are non-pyrogenic.



Lumen	Priming Volume* (mL)	Flow Rate† (mL/hr)
Distal (16 Ga.)	0.51	2100
Medial 1 (14 Ga.)	0.62	3900
Medial 2 (18 Ga.)	0.41	1000
Proximal (18 Ga.)	0.44	1000

* Priming volumes are approximate and are done without accessories.

† Based on 100 cm head height.



LOT LotNumber

YYYY-MM-DD



(01)XXXXXXXXXXXXXX
(17)YYMMDD
(10)LotNumber

CE 0086
Teleflex Medical
IDA Business and
Technology Park
Dublin Road, Athlone
Co. Westmeath, Ireland

Arrow International, Inc.
Subsidiary of Teleflex Incorporated
2400 Bernville Road
Reading, PA 19605 USA

Teleflex®
www.teleflex.com

LBL010390 (08/20)



REF CS-14854-E 8.5 Fr. 4 Lumen 30 cm catheter length .032 inch dia. spring-wire guide

Quad-Lumen Central Venous Catheterization Set with Blue FlexTip® Catheter

