



REF CS-15853-E (IPNXXXXXX)

8.5 Fr. 3 Lumen 20 cm catheter length .032 inch dia. spring-wire guide

Multi-Lumen Central Venous Catheterization Set with Blue FlexTip® Catheter

Contents:

- 1: Multi-Lumen Indwelling Catheter: 8.5 Fr. x 20 cm Radiopaque Polyurethane with Blue FlexTip®, Extension Line Clamps
- 1: Spring-Wire Guide, Marked: .032" (0.81 mm) dia. x 23-5/8" (60 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer
- 1: Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle
- 1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW
- 1: Pressure Transduction Probe
- 1: Arrow® Raulerson Spring-Wire Introduction Syringe: 5 mL
- 1: Tissue Dilator: 9 Fr. (3.0 mm) x 10.2 cm
- 3: Dust Cap: Non-Vented
- 1: SecondSite™ Adjustable Hub: Fastener
- 1: SecondSite™ Adjustable Hub: Catheter Clamp

All components are CE 0086 unless otherwise noted.

Rx only. Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.

California Prop. 65


⚠ **WARNING:** Cancer and Reproductive Harm.

www.P65Warnings.ca.gov

Not made with natural rubber latex.

Fluid path components are non-pyrogenic.

- tw** 帶 Blue FlexTip 導管的多腔中央靜脈插管套裝
- cs** Vicelumenový centrální žilní katetizační set s katetrem Blue FlexTip
- da** Multilumen centralvenekateriseringsset med Blue FlexTip kateter
- nl** Multilumen centraalveneuze katheterisatieset met Blue FlexTip katheter
- fr** Jeu de cathétérisme veineux central à lumières multiples avec cathéter Blue FlexTip
- de** Besteck für viellumigen zentralen Venenkatheter mit Blue FlexTip Katheter
- el** Σετ κεντρικού φλεβικού καθετηριασμού πολλαπλών αυλών με καθετήρα Blue FlexTip
- it** Set per cateterismo della vena centrale multilume con catetere Blue FlexTip
- ja** マルチルーメン中心静脈カテーテルセット-Blue FlexTipカテーテル付
- no** Multilumen sentralt venekateriseringssett med Blue FlexTip kateter
- pl** Zestaw wielokanalowy do cewnikowania żył centralnych z cewnikiem typu Blue FlexTip
- pt** Conjunto de cateterização venosa central com lúmen múltiplo e cateter Blue FlexTip
- ru** Набор с многопросветным центральным венозным катетером с кончиком Blue FlexTip
- sk** Multilúmenová centrálna venózna katetizačná sada s katétrom Blue FlexTip
- sl** Večlumenski centralni venski komplet za katetizacijo s katetrom Blue FlexTip
- es** Conjunto de cateterización venosa central de luz múltiple con catéter Blue FlexTip
- sv** Set för flerkanalig central venkatetrisering med Blue FlexTip-kateter
- tr** Multi-Lümen Santral Venöz Katetizasyon Seti, Blue FlexTip Kateter ile



Lumen	Priming Volume* (ml)	Flow Rate† (ml/hr)
Distal (16 Ga.)	0.47	3530
Medial (14 Ga.)	0.89	7555
Proximal (16 Ga.)	0.47	3290

* Priming volumes are approximate and are done without accessories.

† Based on 100 cm head height.



LOT LotNumber
YYYY-MM-DD



(01)XXXXXXXXXXXXXX
(17)YYMMDD
(10)LotNumber

CE REP
0086
Teleflex Medical
IDA Business and
Technology Park
Dublin Road, Athlone
Co. Westmeath, Ireland

Arrow International, Inc.
Subsidiary of Teleflex Incorporated
2400 Bernville Road
Reading, PA 19605 USA

Teleflex®
www.teleflex.com

LBL010399 (02/19)



REF CS-15853-E 8.5 Fr. 3 Lumen 20 cm catheter length .032 inch dia. spring-wire guide

Multi-Lumen Central Venous Catheterization Set with Blue FlexTip® Catheter

