



**REF** CS-15955 (IPNXXXXXX)

**9.5** Fr. **5** Lumen **20** cm catheter length **.032** inch dia. spring-wire guide

# Multi-Lumen Central Venous Catheterization Set with Blue FlexTip® Catheter

**Contents:**


- 1: Multi-Lumen Indwelling Catheter: 9.5 Fr. x 20 cm Radiopaque Polyurethane with Blue FlexTip®, Extension Line Clamps
- 1: Spring-Wire Guide, Marked: .032" (0.81 mm) dia. x 23-5/8" (60 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer with ECG Mark
- 1: Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle
- 1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW
- 1: Syringe: 5 mL Luer-Slip
- 1: Tissue Dilator: 10 Fr. (3.5 mm) x 10.2 cm
- 5: Dust Cap: Non-Vented
- 1: SecondSite™ Adjustable Hub: Fastener
- 1: SecondSite™ Adjustable Hub: Catheter Clamp

All components are CE 0086 unless otherwise noted.

**Rx only. Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.**

Not made with natural rubber latex.

Fluid path components are non-pyrogenic.



Lumen	Priming Volume* (ml)	Flow Rate† (ml/hr)
Distal (16 Ga.)	0.48	3605
Medial 1 (14 Ga.)	0.57	5390
Medial 2 (18 Ga.)	0.38	1250
Medial 3 (18 Ga.)	0.38	1160
Proximal (18 Ga.)	0.41	1455

\* Priming volumes are approximate and are done without accessories.

† Based on 100 cm head height.

- be** Многопросветный набор для центральной венозной катетеризации с катэтрам Blue FlexTip
- hr** Višelumenski centralni venski kateter Set s Blue FlexTip kateterom
- cs** Vičelumenový centrální žilní katetrizační set s katetrem Blue FlexTip
- da** Multilumen centralvenøs kateteriseringssæt med Blue FlexTip-kateter
- nl** Multilumen centraalveneuze katheterisatieset met Blue FlexTip katheter
- et** Mitme valendikuga tsentraalse veenikateeteri komplekt koos Blue FlexTip kateetriga
- fi** Monilumeninen keskuslaskimokatetrisarja jossa Blue FlexTip-katetri
- fr** Jeu de cathétérisme veineux central à lumières multiples avec cathéter Blue FlexTip
- de** Besteck für viellumigen, zentralen Venenkatheter mit Blue-FlexTip-Katheter
- el** Σετ Καθετηριασμού Κεντρικής φλέβας Πολλαπλούς Αυλούς με Καθετήρα Blue FlexTip
- hu** Többlumenű, centrális vénás katéterezési szett, Blue FlexTip katéterrel
- is** Fjöl lúmena holæðaleggsett með Blue FlexTip æðalegg
- it** Set per cateterismo della vena centrale multilume con catetere Blue FlexTip
- lt** Daugiakanalis centrínės venos kateterizavimo mini rinkinys su „Blue FlexTip“ kateteriu
- no** Multilumen sentralt venekateteriseringssett med Blue FlexTip kateter
- pl** Zestaw wielokanalowy do cewnikowania żył centralnych z cewnikiem typu Blue FlexTip
- pt** Conjunto de Cateterização Venosa Central com Lúmen Múltiplo e Cateter Blue FlexTip
- ro** Set de cateterizare venoasă centrală multi-lumen cu cateter Blue FlexTip
- ru** Набор с многопросветным центральным венозным катетером с кончиком Blue FlexTip
- sr** Višelumenski centralni venski kateterizacioni Set sa Blue FlexTip kateterom
- sk** Multilúmenová centrálna venózna katetrizačná sada s katétrom Blue FlexTip
- sl** Večlumenski centralni venski komplet za katetrizacijo s katetrom Blue FlexTip
- es** Conjunto de cateterización venosa central de luz múltiple con catéter Blue FlexTip
- sv** Set för flerkanalig central venkatetrisering med Blue FlexTip-kateter
- tr** Multi Lümen Santral Venöz Kateterizasyon Seti, Blue FlexTip Kateter ile
- uk** Багатопрсвітний центральний венозний катетер Набір: катетер Blue с кінчиком Blue FlexTip Catheter



**LOT** LotNumber  
YYYY-MM-DD



(01)XXXXXXXXXXXXXX  
(17)YYMMDD  
(10)LotNumber

**CE** REP 0086  
Teleflex Medical  
IDA Business and  
Technology Park  
Dublin Road, Athlone  
Co. Westmeath, Ireland

Arrow International, Inc.  
Subsidiary of Teleflex Incorporated  
2400 Bernville Road  
Reading, PA 19605 USA

**Teleflex**  
www.teleflex.com  
LBL010402 (02/19)



**REF** CS-15955 **9.5** Fr. **5** Lumen **20** cm catheter length **.032** inch dia. spring-wire guide

Multi-Lumen Central Venous Catheterization Set with Blue FlexTip® Catheter

