



REF CS-22853-E (IPNXXXXXX)

8.5 Fr. 3 Lumen 16 cm catheter length .032 inch dia. spring-wire guide

Multi-Lumen Central Venous Catheterization Set with Blue FlexTip® ARROWg+ard Blue® Catheter

Contents:

- 1: Multi-Lumen Indwelling Catheter: 8.5 Fr. x 16 cm Radiopaque Polyurethane with Blue FlexTip®, ARROWg+ard® Antimicrobial Surface Treatment¹, Extension Line Clamps
- 1: Spring-Wire Guide, Marked: .032" (0.81 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End – "J" Tip on Other) with Arrow Advancer
- 1: Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle
- 1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW
- 1: Pressure Transduction Probe
- 1: Arrow® Raulerson Spring-Wire Introduction Syringe: 5 mL
- 1: Tissue Dilator: 9 Fr. (3.0 mm) x 10.2 cm
- 3: Dust Cap: Non-Vented
- 1: SecondSite™ Adjustable Hub: Fastener
- 1: SecondSite™ Adjustable Hub: Catheter Clamp

¹Licensed under US Patent No. 6,706,024.

All components are CE 0086 unless otherwise noted.

Rx only. Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.

California Prop. 65


⚠ **WARNING:** Cancer and Reproductive Harm. www.P65Warnings.ca.gov

Not made with natural rubber latex. Store below 25°C (77°F). Avoid excessive heat above 40°C (104°F).

Fluid path components are non-pyrogenic.

Contraindications: The Arrowg+ard Blue antimicrobial catheter is contraindicated for patients with known hypersensitivity to chlorhexidine, silver sulfadiazine and/or sulfa drugs (refer to product instructions for references).

- cs** Vícelumenový centrální žilní katetizační set s katetrem Blue FlexTip ARROWg+ard Blue
- da** Multilumen centralvenekateteriseringssett med Blue FlexTip ARROWg+ard Blue kateter
- nl** Multi-lumen centraal veneuze katheterisatieset met Blue FlexTip ARROWg+ard Blue kateter
- et** Mitme valendikuga tsentraalse veenikateetri komplekt koos Blue FlexTip ARROWg+ard Blue kateetriga
- fi** Monilumeninen keskuskaslaskimokatetrisarja, jossa sininen FlexTip ARROWg+ard Blue katetri
- fr** Jeu de cathétérisme veineux central à lumières multiples avec cathéter Blue FlexTip ARROWg+ard Blue
- de** Besteck für viellumigen, zentralen Venenkatheter mit Blue FlexTip ARROWg+ard Blue Katheter
- el** Σετ Καθητηριασμού Κεντρικής φλέβας με Πολλαπλούς Αυτούς και Καθετήρα ARROWg+ard Blue με Blue FlexTip
- hu** Többlumenű, centrális vénás katéterezési szett, Blue FlexTip ARROWg+ard Blue katéterrel
- it** Set per cateterismo della vena centrale multilume con catetere Blue FlexTip ARROWg+ard Blue
- no** Multilumen Sentralt venekateteriseringssett med Blue FlexTip ARROWg+ard Blue kateter
- pl** Zestaw do cewnikowania żyły centralnej z cewnikiem wielokanalowym z końcówką Blue FlexTip i powłoką antybakteryjną ARROWg+ard Blue
- pt** Conjunto de Cateterização Venosa Central com Lúmen Múltiplo e Cateter Blue FlexTip, ARROWg+ard Blue
- ru** Набор с многопросветным центральным венозным катетером ARROWg+ard Blue с кончиком Blue FlexTip
- sk** Viacľúmenová centrálna venózna katetizačná súprava s katétrom Blue FlexTip ARROWg+ard Blue
- sl** Večlumenski centralni venski komplet za katetizacijo s katetrom Blue FlexTip ARROWg+ard Blue
- es** Conjunto de cateterización venosa central de luz múltiple con catéter Blue FlexTip ARROWg+ard Blue
- sv** Set för flerkanalig central venkatetrisering med Blue FlexTip ARROWg+ard Blue kateter
- tr** Multi Lümen Santral Venöz Kateterizasyon Seti, Mavi esnek uçlu ARROWg+ard Blue Kateterli



Lumen	Priming Volume* (ml)	Flow Rate† (ml/hr)
Distal (16 Ga.)	0.43	3950
Medial (14 Ga.)	0.78	7630
Proximal (16 Ga.)	0.43	3820

* Priming volumes are approximate and are done without accessories.

† Based on 100 cm head height.



LOT LotNumber
YYYY-MM-DD



(01)XXXXXXXXXXXXXX
(17)YYMMDD
(10)LotNumber



CE 0086
Teleflex Medical
IDA Business and
Technology Park
Dublin Road, Athlone
Co. Westmeath, Ireland

Arrow International, Inc.
Subsidiary of Teleflex Incorporated
2400 Bernville Road
Reading, PA 19605 USA

Teleflex®
www.teleflex.com

LBL010408 (01/19)



REF CS-22853-E 8.5 Fr. 3 Lumen 16 cm catheter length .032 inch dia. spring-wire guide

Multi-Lumen Central Venous Catheterization Set with Blue FlexTip® ARROWg+ard Blue® Catheter

