



REF CS-22855 (IPNXXXXXX)

8.5 Fr. 5 Lumen 16 cm catheter length .032 inch dia. spring-wire guide

Multi-Lumen Central Venous Catheterization Set with Blue FlexTip® ARROWg+ard Blue® Catheter

Contents:

- 1: Multi-Lumen Indwelling Catheter: 8.5 Fr. x 16 cm Radiopaque Polyurethane with Blue FlexTip®, ARROWg+ard® Antimicrobial Surface Treatment¹, Extension Line Clamps
- 1: Spring-Wire Guide, Marked: .032" (0.81 mm) dia. x 23-5/8" (60 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer with ECG Mark
- 1: Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle
- 1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW
- 1: Syringe: 5 mL Luer-Slip
- 1: Tissue Dilator: 9 Fr. (3.0 mm) x 10.2 cm
- 5: Dust Cap: Non-Vented
- 1: SecondSite™ Adjustable Hub: Fastener
- 1: SecondSite™ Adjustable Hub: Catheter Clamp
- 1: Scalpel: #11

¹Licensed under US Patent No. 6,706,024.

All components are CE 0086 unless otherwise noted.


Rx only. Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.

Not made with natural rubber latex. Store below 25°C (77°F). Avoid excessive heat above 40°C (104°F).

Fluid path components are non-pyrogenic.

Contraindications: The Arrowg+ard Blue antimicrobial catheter is contraindicated for patients with known hypersensitivity to chlorhexidine, silver sulfadiazine and/or sulfa drugs (refer to product instructions for references).

- cs** Vícelumenový centrální žilní katetrizační set s katetrem Blue FlexTip ARROWg+ard Blue
- da** Multilumen centralvenekateteriseringssæt med Blue FlexTip ARROWg+ard Blue kateter
- nl** Multi-lumen centraal veneuze katheterisatieset met Blue FlexTip ARROWg+ard Blue katheter
- et** Mitme valendikuga tsentraalse veenikateetri komplekt koos Blue FlexTip ARROWg+ard Blue kateetriga
- fi** Monilumeninen keskuslaskimokatetri sarja, jossa sininen FlexTip ARROWg+ard Blue katetri
- fr** Jeu de cathétérisme veineux central à lumières multiples avec cathéter Blue FlexTip ARROWg+ard Blue
- de** Besteck für viellumigen, zentralen Venenkatheter mit Blue FlexTip ARROWg+ard Blue Katheter
- el** Σετ Καθετηριασμού Κεντρικής φλέβας με Πολλαπλούς Αυτούς και Καθετήρα ARROWg+ard Blue με Blue FlexTip
- hu** Többlumenű, centrális vénás katéterezési szett, Blue FlexTip ARROWg+ard Blue katéterrel
- it** Set per cateterismo della vena centrale multilume con catetere Blue FlexTip ARROWg+ard Blue
- no** Multilumen Sentralt venekateteriseringssæt med Blue FlexTip ARROWg+ard Blue kateter
- pl** Zestaw do cewnikowania żyły centralnej z cewnikiem wielokanalowym z końcówką Blue FlexTip i powłoką antybakteryjną ARROWg+ard Blue
- pt** Conjunto de Cateterização Venosa Central com Lúmen Múltiplo e Cateter Blue FlexTip, ARROWg+ard Blue
- sk** Viacelumenová centrálna venózna katetrizačná súprava s katetrom Blue FlexTip ARROWg+ard Blue
- es** Conjunto de cateterización venosa central de luz múltiple con catéter Blue FlexTip ARROWg+ard Blue
- sv** Set för flerkanalig central venkatetrering med Blue FlexTip ARROWg+ard Blue kateter



Lumen	Priming Volume* (ml)	Flow Rate† (ml/hr)
Distal (16 Ga.)	0.58	3650
Medial 1 (14 Ga.)	0.71	6590
Medial 2 (18 Ga.)	0.36	1260
Medial 3 (18 Ga.)	0.38	1430
Proximal (18 Ga.)	0.4	1460

* Priming volumes are approximate and are done without accessories.

† Based on 100 cm head height.



LOT LotNumber
YYYY-MM-DD



(01)XXXXXXXXXXXXXX
(17)YYMMDD
(10)LotNumber

CE REP Teleflex Medical
IDA Business and
Technology Park
Dublin Road, Athlone
Co. Westmeath, Ireland
0086

Arrow International, Inc.
Subsidiary of Teleflex Incorporated
2400 Bernville Road
Reading, PA 19605 USA

Teleflex®
www.teleflex.com

LBL010410 (02/19)



REF CS-22855 8.5 Fr. 5 Lumen 16 cm catheter length .032 inch dia. spring-wire guide

Multi-Lumen Central Venous Catheterization Set with Blue FlexTip® ARROWg+ard Blue® Catheter

