

**ARROW**



**REF CS-25402**

(IPNXXXXXX)

**4** Fr. | **2** Lumen | **8** cm catheter length | **.018** inch dia. spring-wire guide

# Pediatric Two-Lumen Central Venous Catheterization Set with Blue FlexTip® ARROWg+ard Blue® Catheter

**Contents:**

- 1: Two-Lumen Indwelling Catheter: 4 Fr. x 8 cm Radiopaque Polyurethane with Blue FlexTip®, ARROWg+ard® Antimicrobial Surface Treatment<sup>1</sup>, Extension Line Clamps
- 1: Spring-Wire Guide, Marked: .018" (0.46 mm) dia. x 17-3/4"(45 cm) (Straight Soft Tip on One End – "J" Tip on Other) with Arrow Advancer
- 1: Catheter: 22 Ga. x 1-3/4" (4.45 cm) Radiopaque over 25 Ga. RW Introducer Needle
- 1: Introducer Needle: 21 Ga. x 1-1/2" (3.81 cm) TW
- 1: Syringe: 5 mL Luer-Slip
- 1: Tissue Dilator: 5 Fr. (1.8 mm) x 7.5 cm
- 1: Tissue Dilator: 5 Fr. (1.8 mm) x 4.1 cm
- 2: Dust Cap: Non-Vented
- 1: SecondSite™ Adjustable Hub: Catheter Clamp
- 1: SecondSite™ Adjustable Hub: Fastener

<sup>1</sup>Licensed under US Patent No. 6,706,024.

All components are CE 0086 unless otherwise noted.

**Rx only. Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.**

California Prop. 65


⚠ **WARNING:** Cancer and Reproductive Harm. [www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov)

Not made with natural rubber latex. Store below 25°C (77°F). Avoid excessive heat above 40°C (104°F).

Fluid path components are non-pyrogenic.

**Contraindications: The ARROWg+ard Blue® antimicrobial catheter is contraindicated for patients with known hypersensitivity to chlorhexidine, silver sulfadiazine and/or sulfa drugs (refer to product instructions for references).**

- fr** Jeu de cathétérisme veineux central pédiatrique à deux lumières avec cathéter Blue FlexTip ARROWg+ard Blue
- de** Besteck für zweilumigen, zentralen Venenkatheter für Kinder mit Blue FlexTip ARROWg+ard Blue Katheter
- it** Set per cateterismo della vena centrale a due lumi pediatrico con catetere Blue FlexTip ARROWg+ard Blue
- pl** Pediatriczny zestaw do cewnikowania żył centralnych, z 2-kanalowym cewnikiem z końcówką Blue FlexTip i powłoką antybakteryjną ARROWg+ard Blue
- pt** Conjunto Pediátrico de Cateterização Venosa Central com Lúmen Duplo e Cateter Blue FlexTip, ARROWg+ard Blue
- ru** Детский набор для центральной венозной катетеризации с двухпросветным катетером ARROWg+ard и наконечником Blue FlexTip
- sl** Dvolumenski centralni venski komplet za katetrijazijo s katetrom Blue FlexTip ARROWg+ard Blue za uporabo v pediatrji
- es** Conjunto de cateterización venosa central de dos luces para uso pediátrico con catéter Blue FlexTip ARROWg+ard Blue
- sv** Set för pediatriisk tvåkanalig central venkatetrisering med Blue FlexTip ARROWg+ard Blue kateter
- tr** Blue FlexTip ARROWg+ard Blue Kateter İçeren Pediyatrik İki Lümenli Santral Venöz Katetrijasyon Seti



Lumen	Priming Volume* (ml)	Flow Rate† (ml/hr)
Distal (22 Ga.)	0.3	880
Proximal (22 Ga.)	0.30	480

\* Priming volumes are approximate and are done without accessories.

† Based on 100 cm head height.



**LOT** LotNumber  
YYYY-MM-DD



(01)XXXXXXXXXXXXXX  
(17)YYMMDD  
(10)LotNumber

**CE** 0086  
Teleflex Medical  
IDA Business and  
Technology Park  
Dublin Road, Athlone  
Co. Westmeath, Ireland

Arrow International, Inc.  
Subsidiary of Teleflex Incorporated  
2400 Bernville Road  
Reading, PA 19605 USA

**Teleflex**  
www.teleflex.com

LBL010412 (11/18)



**REF CS-25402**

**4** Fr. | **2** Lumen | **8** cm catheter length | **.018** inch dia. spring-wire guide

**Pediatric Two-Lumen Central Venous Catheterization Set with Blue FlexTip® ARROWg+ard Blue® Catheter**

