



REF CS-25855-E (IPNXXXXXX)

8.5 Fr. 5 Lumen 20 cm catheter length .032 inch dia. spring-wire guide

# Multi-Lumen Central Venous Catheterization Set with Blue FlexTip® ARROWg+ard Blue® Catheter

**Contents:**

- 1: Multi-Lumen Indwelling Catheter: 8.5 Fr. x 20 cm Radiopaque Polyurethane with Blue FlexTip®, ARROWg+ard® Antimicrobial Surface Treatment<sup>1</sup>, Extension Line Clamps
- 1: Spring-Wire Guide, Marked: .032" (0.81 mm) dia. x 23-5/8" (60 cm) (Straight Soft Tip on One End – "J" Tip on Other) with Arrow Advancer with ECG Mark
- 1: Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle
- 1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW
- 1: Pressure Transduction Probe
- 1: Arrow® Raulerson Spring-Wire Introduction Syringe: 5 mL
- 1: Tissue Dilator: 9 Fr. (3.0 mm) x 10.2 cm
- 5: Dust Cap: Non-Vented
- 1: SecondSite™ Adjustable Hub: Fastener
- 1: SecondSite™ Adjustable Hub: Catheter Clamp

<sup>1</sup>Licensed under US Patent No. 6,706,024.

All components are CE 0086 unless otherwise noted.


**Rx only. Warning:** Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.

Not made with natural rubber latex. Store below 25°C (77°F). Avoid excessive heat above 40°C (104°F).

Fluid path components are non-pyrogenic.

**Contraindications:** The Arrowg+ard Blue antimicrobial catheter is contraindicated for patients with known hypersensitivity to chlorhexidine, silver sulfadiazine and/or sulfa drugs (refer to product instructions for references).

- cs** Vícelumenový centrální žilní katetrizační set s katetrem Blue FlexTip ARROWg+ard Blue
- da** Multilumen centralvenekateriseringsset med Blue FlexTip ARROWg+ard Blue kateter
- nl** Multi-lumen centraal veneuze katheterisatieset met Blue FlexTip ARROWg+ard Blue katheter
- et** Mitme valendikuga tsentraalse veenikateetri komplekt koos Blue FlexTip ARROWg+ard Blue kateetriga
- fi** Monilumeninen keskuslaskimokatetri sarja, jossa sininen FlexTip ARROWg+ard Blue katetri
- fr** Jeu de cathétérisme veineux central à lumières multiples avec cathéter Blue FlexTip ARROWg+ard Blue
- de** Besteck für viellumigen, zentralen Venenkatheter mit Blue FlexTip ARROWg+ard Blue Katheter
- el** Σετ Καθετηριασμού Κεντρικής φλέβας με Πολλαπλούς Αυτούς και Καθετήρα ARROWg+ard Blue με Blue FlexTip
- hu** Többlumenű, centrális vénás katéterezési szett, Blue FlexTip ARROWg+ard Blue katéterrel
- it** Set per cateterismo della vena centrale multilume con catetere Blue FlexTip ARROWg+ard Blue
- no** Multilumen Sentralt venekateriseringssett med Blue FlexTip ARROWg+ard Blue kateter
- pl** Zestaw do cewnikowania żyły centralnej z cewnikiem wielokanalowym z końcówką Blue FlexTip i powłoką antybakteryjną ARROWg+ard Blue
- pt** Conjunto de Cateterização Venosa Central com Lúmen Múltiplo e Cateter Blue FlexTip, ARROWg+ard Blue
- sk** Viacelumenová centrálna venózna katetrizačná súprava s katétrom Blue FlexTip ARROWg+ard Blue
- es** Conjunto de cateterización venosa central de luz múltiple con catéter Blue FlexTip ARROWg+ard Blue
- sv** Set för flerkanalig central venkatetrering med Blue FlexTip ARROWg+ard Blue kateter



Lumen	Priming Volume* (ml)	Flow Rate† (ml/hr)
Distal (16 Ga.)	0.63	3240
Medial 1 (14 Ga.)	0.78	5720
Medial 2 (18 Ga.)	0.38	1070
Medial 3 (18 Ga.)	0.4	1220
Proximal (18 Ga.)	0.42	1190

\* Priming volumes are approximate and are done without accessories.

† Based on 100 cm head height.



LOT LotNumber  
YYYY-MM-DD



(01)XXXXXXXXXXXXXX  
(17)YYMMDD  
(10)LotNumber



CE 0086  
Teleflex Medical  
IDA Business and  
Technology Park  
Dublin Road, Athlone  
Co. Westmeath, Ireland

Arrow International, Inc.  
Subsidiary of Teleflex Incorporated  
2400 Bernville Road  
Reading, PA 19605 USA

Teleflex®  
www.teleflex.com

LBL010422 (01/19)



REF CS-25855-E 8.5 Fr. 5 Lumen 20 cm catheter length .032 inch dia. spring-wire guide

Multi-Lumen Central Venous Catheterization Set with Blue FlexTip® ARROWg+ard Blue® Catheter

