



REF CS-26553 (IPNXXXXXX)

5.5 Fr. 3 Lumen 13 cm catheter length .018 inch dia. spring-wire guide

Pediatric Multi-Lumen Central Venous Catheterization Set with Blue FlexTip® ARROWg+ard Blue® Catheter

Contents:

- 1: Multi-Lumen Indwelling Catheter: 5.5 Fr. x 13 cm Radiopaque Polyurethane with Blue FlexTip®, ARROWg+ard® Antimicrobial Surface Treatment†, Extension Line Clamps
- 1: Spring-Wire Guide, Marked: .018" (0.46 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer
- 1: Catheter: 22 Ga. x 1-3/4" (4.45 cm) Radiopaque over 25 Ga. RW Introducer Needle
- 1: Introducer Needle: 21 Ga. x 1-1/2" (3.81 cm) TW
- 1: Syringe: 5 mL Luer-Slip
- 1: Tissue Dilator: 6 Fr. (2.0 mm) x 7.6 cm
- 1: Tissue Dilator: 6 Fr. (2.0 mm) x 4.1 cm
- 3: Dust Cap: Non-Vented
- 1: SecondSite™ Adjustable Hub: Fastener
- 1: SecondSite™ Adjustable Hub: Catheter Clamp

†Licensed under US Patent No. 6,706,024.

All components are CE 0086 unless otherwise noted.

Rx only. Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.

California Prop. 65


⚠ **WARNING:** Cancer and Reproductive Harm. www.P65Warnings.ca.gov

Not made with natural rubber latex. Store below 25°C (77°F). Avoid excessive heat above 40°C (104°F).

Fluid path components are non-pyrogenic.

Contraindications: The Arrowg+ard Blue antimicrobial catheter is contraindicated for patients with known hypersensitivity to chlorhexidine, silver sulfadiazine and/or sulfa drugs (refer to product instructions for references).

- fr** Jeu de cathétérisme veineux central pédiatrique à lumières multiples avec cathéter Blue FlexTip ARROWg+ard Blue
- de** Besteck für viellumigen, zentralen Venenkatheter für Kinder mit Blue-FlexTip-ARROWg+ard-Blue-Katheter
- it** Set per cateterismo della vena centrale multilume pediatrico con catetere Blue FlexTip ARROWg+ard Blue
- pl** Wielokanałowy zestaw pediatryczny do cewnikowania żył centralnych z cewnikami Blue FlexTip ARROWg+ard Blue
- pt** Conjunto Pediátrico de Cateterização Venosa Central com Lúmen Múltiplo e Cateter Blue FlexTip, ARROWg+ard Blue
- ru** Набор для многопросветной катетеризации центральных вен у детей с катетером Blue FlexTip ARROWg+ard Blue
- sl** Pediatrični komplet za centralno venosko katetizacijo z več svetlinami s katetrom Blue FlexTip ARROWg+ard Blue
- es** Conjunto de cateterización venosa central de luz múltiple para uso pediátrico con catéter Blue FlexTip ARROWg+ard Blue
- sv** Set för pediatriisk flerkanalig central venkatetrisering med Blue FlexTip ARROWg+ard Blue-kateter
- tr** Pediyatrik Multi-Lümen Santral Venöz Kateterizasyon Seti, Blue FlexTip ARROWg+ard Blue Kateter ile



Lumen	Priming Volume* (ml)	Flow Rate† (ml/hr)
Distal (20 Ga.)	0.30	1220
Medial (22 Ga.)	0.3	560
Proximal (22 Ga.)	0.4	670

* Priming volumes are approximate and are done without accessories.

† Based on 100 cm head height.



LOT LotNumber
YYYY-MM-DD



(01)XXXXXXXXXXXXXX
(17)YYMMDD
(10)LotNumber



CE 0086
Teleflex Medical
IDA Business and
Technology Park
Dublin Road, Athlone
Co. Westmeath, Ireland

Arrow International, Inc.
Subsidiary of Teleflex Incorporated
2400 Bernville Road
Reading, PA 19605 USA

Teleflex®
www.teleflex.com

LBL010424 (02/19)



REF CS-26553 5.5 Fr. 3 Lumen 13 cm catheter length .018 inch dia. spring-wire guide

Pediatric Multi-Lumen Central Venous Catheterization Set with Blue FlexTip® ARROWg+ard Blue® Catheter

