



REF CS-42703-E (IPNXXXXXX)

7 Fr. 3 Lumen 16 cm catheter length .032 inch dia. spring-wire guide

# Arrowg+ard Blue PLUS® Multi-Lumen CVC Set

**Contents:**

- 1: Arrowg+ard Blue Plus®<sup>1</sup> Three-Lumen Catheter: 7 Fr. (2.5 mm OD) x 16 cm
- 1: Spring-Wire Guide, Marked: .032" (0.81 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer
- 1: Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle
- 1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW
- 1: Pressure Transduction Probe
- 1: Arrow® Raulerson Spring-Wire Introduction Syringe: 5 mL
- 1: Tissue Dilator: 8.5 Fr. (2.8 mm) x 10.2 cm
- 3: Dust Cap: Non-Vented
- 1: SecondSite™ Adjustable Hub: Fastener
- 1: SecondSite™ Adjustable Hub: Catheter Clamp

<sup>1</sup>Licensed under US Patent No. 7,329,412.

All components are CE 0086 unless otherwise noted.

**Rx only. Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.**

Not made with natural rubber latex. Store below 25°C (77°F). Avoid excessive heat above 40°C (104°F).

Fluid path components are non-pyrogenic.

**Contraindications: The Arrowg+ard Blue Plus antimicrobial catheter is contraindicated for patients with known hypersensitivity to chlorhexidine, silver sulfadiazine and/or sulfa drugs.**



Lumen	Priming Volume* (mL)	Flow Rate† (mL/hr)
Distal (16 Ga.)	0.38	3400
Medial (18 Ga.)	0.37	1800
Proximal (18 Ga.)	0.39	1900

\* Priming volumes are approximate and are done without accessories.

† Based on 100 cm head height.

- cs** **Vícelenový CVC set Arrowg+ard Blue PLUS**  
Kontraindikace: Použití antimikrobiálního katetru Arrowg+ard Blue Plus je kontraindikováno u pacientů se známou hypersenzitivitou na chlorhexidin, sulfadiazin stříbra a/nebo sulfonamidy.
- nl** **Arrowg+ard Blue PLUS multilumen centraal-veneuze katheterset** Contra-indicaties: De Arrowg+ard Blue Plus antimicrobiële katheter is gecontra-indiceerd voor patiënten met bekende overgevoeligheid voor chloorhexidine, zilver-sulfadiazine en/of sulfonamidebevattende geneesmiddelen.
- fr** **Jeu de cathéters veineux centraux à lumières multiples Arrowg+ard Blue PLUS** Contre-indications : Le cathéter antimicrobien Arrowg+ard Blue Plus est contre-indiqué pour les patients qui présentent une hypersensibilité connue à la chlorhexidine, à la sulfadiazine d'argent et/ou aux sulfamides.
- de** **Set mit mehrlumigen Arrowg+ard Blue PLUS ZVK** Kontraindikationen: Der Arrowg+ard Blue Plus antimikrobielle Katheter ist bei Patienten mit einer bekannten Überempfindlichkeit gegenüber Chlorhexidin, Silber-Sulfadiazin und/oder Sulfa-Medikamenten kontraindiziert.
- it** **Set CVC multilume Arrowg+ard Blue PLUS** Controindicazioni – Il catetere antimicrobico Arrowg+ard Blue Plus è controindicato nei pazienti con ipersensibilità nota alla clorexidina, alla sulfadiazina argenticca e/o ai sulfamidici.
- pt** **Conjunto de CVC de Lúmen Múltiplo Arrowg+ard Blue PLUS** Contraindicações: O cateter antimicrobiano Arrowg+ard Blue Plus está contraindicado em doentes com hipersensibilidade conhecida à clorexidina, à sulfadiazina de prata e/ou a fármacos sulfa.
- es** **Conjunto de CVC de luz múltiple Arrowg+ard Blue PLUS** Contraindicaciones: El catéter antimicrobiano Arrowg+ard Blue Plus está contraindicado en pacientes con hipersensibilidad conocida a la clorhexidina, la sulfadiazina de plata o las sulfamidas.
- sv** **Arrowg+ard Blue PLUS flerlumen CVK-set** Kontraindikationer: Arrowg+ard Blue Plus antimikrobiell kateter är kontraindicerad för patienter med känd överkänslighet mot klorhexidin, silversulfadiazin och/eller sulfaläkemedel.



Arrow International, Inc.  
Subsidiary of Teleflex Incorporated  
2400 Bernville Road  
Reading, PA 19605 USA

LOT LotNumber

YYYY-MM-DD



(01)XXXXXXXXXXXXXX  
(17)YYMMDD  
(10)LotNumber

CE 0086  
Teleflex Medical  
IDA Business and  
Technology Park  
Dublin Road, Athlone  
Co. Westmeath, Ireland

**Teleflex®**  
www.teleflex.com

LBL010426 (08/20)



REF CS-42703-E

7 Fr. 3 Lumen 16 cm catheter length .032 inch dia. spring-wire guide

# Arrowg+ard Blue PLUS® Multi-Lumen CVC Set

