



REF CS-42854-E (IPNXXXXXX)

8.5^{Fr.} 4^{Lumen} 16^{cm} catheter length .032^{inch dia.} spring-wire guide

ARROWg+ard Blue PLUS® Quad-Lumen CVC Set

Contents:

- 1: Quad-Lumen Indwelling Catheter: 8.5 Fr. x 16 cm Radiopaque Polyurethane with Blue FlexTip®, ARROWg+ard Blue Plus® Antimicrobial Surface Treatment[†], Extension Line Clamps
- 1: Spring-Wire Guide, Marked: .032" (0.81 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer
- 1: Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle
- 1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW
- 1: Pressure Transduction Probe
- 1: Arrow® Raulerson Spring-Wire Introduction Syringe: 5 mL
- 1: Tissue Dilator: 9 Fr. (3.0 mm) x 10.2 cm
- 4: Dust Cap: Non-Vented
- 1: SecondSite™ Adjustable Hub: Fastener
- 1: SecondSite™ Adjustable Hub: Catheter Clamp

[†]Licensed under US Patent Nos. 6,706,024; 6,872,195 and 7,329,412.

All components are CE 0086 unless otherwise noted.


- cs** Čtyřlumenový CVC set ARROWg+ard Blue PLUS
- nl** ARROWg+ard Blue PLUS centraal-veneuze katheterset met vier lumina
- fr** Jeu de cathéters veineux centraux à quatre lumières ARROWg+ard Blue PLUS
- de** Set mit vierlumigen ARROWg+ard Blue PLUS ZVK
- it** Set CVC quadrilume ARROWg+ard Blue PLUS
- pt** Conjunto de CVC de Lúmen quádruplo ARROWg+ard Blue PLUS
- es** Conjunto de CVC de luz quádruple ARROWg+ard Blue PLUS
- sv** ARROWg+ard Blue PLUS CVK-set med fyra lumen

Rx only. Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.

Not made with natural rubber latex. Store below 25°C (77°F). Avoid excessive heat above 40°C (104°F).

Fluid path components are non-pyrogenic.

Contraindications: The Arrowg+ard Blue PLUS antimicrobial catheter is contraindicated for patients with known hypersensitivity to chlorhexidine, silver sulfadiazine and/or sulfa drugs (refer to product instructions for references).



Lumen	Priming Volume* (ml)	Flow Rate† (ml/hr)
Distal (16 Ga.)	0.39	3200
Medial 1 (14 Ga.)	0.47	5400
Medial 2 (18 Ga.)	0.36	1600
Proximal (18 Ga.)	0.38	1800

* Priming volumes are approximate and are done without accessories.

† Based on 100 cm head height.



LOT LotNumber
YYYY-MM-DD



(01)XXXXXXXXXXXXXX
(17)YYMMDD
(10)LotNumber



EC REP Teleflex Medical
IDA Business and
Technology Park
Dublin Road, Athlone
Co. Westmeath, Ireland



Arrow International, Inc.
Subsidiary of Teleflex Incorporated
2400 Bernville Road
Reading, PA 19605 USA

Teleflex®
www.teleflex.com

LBL010430 (01/19)



REF CS-42854-E

8.5^{Fr.}

4^{Lumen}

16^{cm} catheter length

.032^{inch dia.} spring-wire guide

ARROWg+ard Blue PLUS® Quad-Lumen CVC Set

