

ARROW



REF CV-15955 (IPNXXXXXX)

9.5 Fr. | **5** Lumen | **20** cm catheter length | **.032** inch dia. spring-wire guide

Multi-Lumen Central Venous Catheterization Set with Blue FlexTip® Catheter

Contents:

- 1: Multi-Lumen Indwelling Catheter: 9.5 Fr. x 20 cm Radiopaque Polyurethane with Blue FlexTip®, Extension Line Clamps
- 1: Spring-Wire Guide, Marked: .032" (0.81 mm) dia. x 23-5/8" (60 cm) (Straight Soft Tip on One End – "J" Tip on Other) with Arrow Advancer with ECG Mark
- 1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW
- 1: Syringe: 5 mL Luer-Slip
- 1: Tissue Dilator: 10 Fr. (3.5 mm) x 10.2 cm
- 5: Dust Cap: Non-Vented
- 1: SecondSite™ Adjustable Hub: Fastener
- 1: SecondSite™ Adjustable Hub: Catheter Clamp


All components are CE 0086 unless otherwise noted.

Rx only. Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.

Not made with natural rubber latex.

Fluid path components are non-pyrogenic.

- be** Многопросветный набор для центральной венозной катетеризации с катэтрам Blue FlexTip
- hr** Višelumenski centralni venski kateter Set s Blue FlexTip kateterom
- cs** Višelumenový centrální žilní katetrizační set s katetrem Blue FlexTip
- da** Multilumen centralvenøs kateteriseringssæt med Blue FlexTip-kateter
- nl** Multilumen centraalveneuze katheterisatieset met Blue FlexTip katheter
- et** Mitme valendikuga tsentraalse veenikateeteri komplekt koos Blue FlexTip kateetriga
- fi** Monilumeninen keskuslaskimokatetrisarja jossa Blue FlexTip-katetri
- fr** Jeu de cathétérisme veineux central à lumières multiples avec cathéter Blue FlexTip
- de** Besteck für viellumigen, zentralen Venenkatheter mit Blue-FlexTip-Katheter
- el** Σετ Καθετηριασμού Κεντρικής φλέβας Πολλαπλούς Αιολούς με Καθετήρα Blue FlexTip
- hu** Többlumenű, centrális vénás katéterezési szett, Blue FlexTip katéterrel
- is** Fjöllumena holæðaleggsett með Blue FlexTip æðalegg
- it** Set per cateterismo della vena centrale multilume con catetere Blue FlexTip
- lt** Daugiakanalės centrinės venos kateterizavimo mini rinkinys su „Blue FlexTip“ kateteriu
- no** Multilumen sentralt venekateteriseringssett med Blue FlexTip kateter
- pl** Zestaw wielokanalowy do cewnikowania żył centralnych z cewnikiem typu Blue FlexTip
- pt** Conjunto de Cateterização Venosa Central com Lúmen Múltiplo e Cateter Blue FlexTip
- ro** Set de cateterizare venoasă centrală multi-lumen cu cateter Blue FlexTip
- ru** Набор с многопросветным центральным венозным катетером с кончиком Blue FlexTip
- sr** Višelumenski centralni venski kateterizacioni Set sa Blue FlexTip kateterom
- sk** Multilúmenová centrálna venózna katetrizačná sada s katétrom Blue FlexTip
- sl** Večlumenski centralni venski komplet za katetrizacijo s katetrom Blue FlexTip
- es** Conjunto de cateterización venosa central de luz múltiple con catéter Blue FlexTip
- sv** Set för flerkanalig central venkatetrisering med Blue FlexTip-kateter
- tr** Multi Lümen Santral Venöz Kateterizasyon Seti, Blue FlexTip Kateter ile
- uk** Багатопрсвітний центральний венозний катетер Набір: катетер Blue с кінчиком Blue FlexTip Catheter



Lumen	Priming Volume* (ml)	Flow Rate† (ml/hr)
Distal (16 Ga.)	0.48	3605
Medial 1 (14 Ga.)	0.57	5390
Medial 2 (18 Ga.)	0.38	1250
Medial 3 (18 Ga.)	0.38	1160
Proximal (18 Ga.)	0.41	1455

* Priming volumes are approximate and are done without accessories.

† Based on 100 cm head height.



LOT LotNumber
YYYY-MM-DD



(01)XXXXXXXXXXXXXX
(17)YYMMDD
(10)LotNumber

CE REP Teleflex Medical
IDA Business and
Technology Park
Dublin Road, Athlone
Co. Westmeath, Ireland
0086

Arrow International, Inc.
Subsidiary of Teleflex Incorporated
2400 Bernville Road
Reading, PA 19605 USA

Teleflex
www.teleflex.com
LBL010447 (11/18)



REF CV-15955

9.5 Fr. | **5** Lumen | **20** cm catheter length | **.032** inch dia. spring-wire guide

Multi-Lumen Central Venous Catheterization Set with Blue FlexTip® Catheter

