



REF CV-42703-E (IPNXXXXXX)

7 Fr. 3 Lumen 16 cm catheter length .032 inch dia. spring-wire guide

Multi-Lumen Central Venous Catheterization Set with ARROWg+ard Blue PLUS® Catheter

Contents:

- 1: Multi-Lumen Indwelling Catheter: 7 Fr. x 16 cm Radiopaque Polyurethane with Blue FlexTip®, ARROWg+ard Blue Plus® Antimicrobial Surface Treatment¹, Extension Line Clamps
- 1: Spring-Wire Guide, Marked: .032" (0.81 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer
- 1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW
- 1: Pressure Transduction Probe
- 1: Arrow® Raulerson Spring-Wire Introduction Syringe: 5 mL
- 1: Tissue Dilator: 8.5 Fr. (2.8 mm) x 10.2 cm
- 3: Dust Cap: Non-Vented
- 1: SecondSite™ Adjustable Hub: Fastener
- 1: SecondSite™ Adjustable Hub: Catheter Clamp

¹Licensed under US Patent Nos. 6,706,024, 6,872,195 and 7,329,412.

All components are CE 0086 unless otherwise noted.


Rx only. Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.

Not made with natural rubber latex. Store below 25°C (77°F). Avoid excessive heat above 40°C (104°F).

Fluid path components are non-pyrogenic.

Contraindications: The Arrowg+ard Blue PLUS antimicrobial catheter is contraindicated for patients with known hypersensitivity to chlorhexidine, silver sulfadiazine and/or sulfa drugs (refer to product instructions for references).

- cs** Vicolumenový CVC set ARROWg+ard Blue PLUS
- nl** ARROWg+ard Blue PLUS multilumen centraal-veneuze katheterset
- fr** Jeu de cathéters veineux centraux à lumières multiples ARROWg+ard Blue PLUS
- de** Set mit mehrlumigen ARROWg+ard Blue PLUS ZVK
- it** Set CVC multilume ARROWg+ard Blue PLUS
- pl** Wielokanałowy zestaw do cewnikowania żył centralnych z cewnikiem ARROWg+ard Blue PLUS
- pt** Conjunto de CVC de Lúmen Múltiplo ARROWg+ard Blue PLUS
- ru** Комплект для катетеризации главных вен с многопросветным катетером ARROWg+ard Blue PLUS
- sl** Komplet za centralno venosko katetризacijo z več svetlinami s katetrom ARROWg+ard Blue PLUS
- es** Conjunto de CVC de luz múltiple ARROWg+ard Blue PLUS
- sv** ARROWg+ard Blue PLUS flerlumen CVK-set
- tr** Multi-Lümen Santral Venöz Kateterezasyon Seti, ARROWg+ard Blue PLUS Kateter ile



Lumen	Priming Volume* (mL)	Flow Rate† (mL/hr)
Distal (16 Ga.)	0.39	3400
Medial (18 Ga.)	0.35	1800
Proximal (18 Ga.)	0.37	1900

* Priming volumes are approximate and are done without accessories.

† Based on 100 cm head height.



LOT LotNumber
YYYY-MM-DD



(01)XXXXXXXXXXXXXX
(17)YYMMDD
(10)LotNumber

CE REP 0086
Teleflex Medical
IDA Business and
Technology Park
Dublin Road, Athlone
Co. Westmeath, Ireland

Arrow International, Inc.
Subsidiary of Teleflex Incorporated
2400 Bernville Road
Reading, PA 19605 USA

Teleflex®
www.teleflex.com

LBL010454 (10/19)



REF CV-42703-E

7 Fr. 3 Lumen 16 cm catheter length .032 inch dia. spring-wire guide

Multi-Lumen Central Venous Catheterization Set with ARROWg+ard Blue PLUS® Catheter

