



REF CV-45854 (IPNXXXXXX)

8.5 Fr. 4 Lumen 20 cm catheter length .032 inch dia. spring-wire guide

# Arrowg+ard Blue Plus® Four-Lumen CVC

**Contents:**

- 1: Arrowg+ard Blue Plus®<sup>1</sup> Four-Lumen Catheter: 8.5 Fr. (3.0 mm OD) x 20 cm
- 1: Spring-Wire Guide, Marked: .032" (0.81 mm) dia. x 23-5/8" (60 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer with ECG Mark
- 1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW
- 1: Syringe: 5 mL Luer-Slip
- 1: Tissue Dilator: 9 Fr. (3.0 mm) x 10.2 cm
- 4: Dust Cap: Non-Vented
- 1: SecondSite™ Adjustable Hub: Fastener
- 1: SecondSite™ Adjustable Hub: Catheter Clamp

<sup>1</sup>Licensed under US Patent No. 7,329,412.

All components are CE 0086 unless otherwise noted.

**Rx only. Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.**

Not made with natural rubber latex. Store below 25°C (77°F). Avoid excessive heat above 40°C (104°F).

Fluid path components are non-pyrogenic.

**Contraindications: The Arrowg+ard Blue Plus antimicrobial catheter is contraindicated for patients with known hypersensitivity to chlorhexidine, silver sulfadiazine and/or sulfa drugs.**

**cs Čtyřlumenový CVC Arrowg+ard Blue Plus** Kontraindikace:

Použití antimikrobiálního katetru Arrowg+ard Blue Plus je kontraindikováno u pacientů se známou hypersenzitivitou na chlorhexidin, sulfadiazin stříbra a/nebo sulfonamidy.

**nl Arrowg+ard Blue Plus-CVC met vier lumina** Contra-

indicaties: De Arrowg+ard Blue Plus antimicrobiële katheter is gecontra-indiceerd voor patiënten met bekende overgevoeligheid voor chlorhexidine, zilver-sulfadiazine en/of sulfonamidebevattende geneesmiddelen.

**fr CVC quadruple lumière Arrowg+ard Blue Plus** Contre-

indications : Le cathéter antimicrobien Arrowg+ard Blue Plus est contre-indiqué pour les patients qui présentent une hypersensibilité connue à la chlorhexidine, à la sulfadiazine d'argent et/ou aux sulfamides.

**de Arrowg+ard Blue Plus vierlumiger ZVK** Kontraindikationen:

Der Arrowg+ard Blue Plus antimikrobielle Katheter ist bei Patienten mit einer bekannten Überempfindlichkeit gegenüber Chlorhexidin, Silber-Sulfadiazin und/oder Sulfa-Medikamenten kontraindiziert.

**it CVC a quattro lumi Arrowg+ard Blue Plus**

Controindicazioni – Il catetere antimicrobico Arrowg+ard Blue Plus è controindicato nei pazienti con ipersensibilità nota alla clorexidina, alla sulfadiazina argintica e/o ai sulfamidici.

**pt CVC de quatro lúmenes Arrowg+ard Blue Plus**

Contraíndicações: O cateter antimicrobiano Arrowg+ard Blue Plus está contraíndicado em doentes com hipersensibilidade conhecida à clorexidina, à sulfadiazina de prata e/ou a fármacos sulfa.

**es CVC de cuatro luces Arrowg+ard Blue Plus**

Contraíndicaciones: El catéter antimicrobiano Arrowg+ard Blue Plus está contraíndicado en pacientes con hipersensibilidad conocida a la clorhexidina, la sulfadiazina de plata o las sulfamidas.

**sv Arrowg+ard Blue Plus fyrlumen-CVK** Kontraindikationer:

Arrowg+ard Blue Plus antimikrobiell kateter är kontraindicerad för patienter med känd överkänslighet mot klorhexidin, silversulfadiazin och/eller sulfaläkemedel.



Lumen	Priming Volume* (mL)	Flow Rate† (mL/hr)
Distal (16 Ga.)	0.40	2700
Medial 1 (14 Ga.)	0.50	5100
Medial 2 (18 Ga.)	0.36	1400
Proximal (18 Ga.)	0.39	1600

\* Priming volumes are approximate and are done without accessories.

† Based on 100 cm head height.



LOT LotNumber

YYYY-MM-DD



(01)XXXXXXXXXXXXXX  
(17)YYMMDD  
(10)LotNumber

CE 0086  
Teleflex Medical  
IDA Business and  
Technology Park  
Dublin Road, Athlone  
Co. Westmeath, Ireland

Arrow International, Inc.  
Subsidiary of Teleflex Incorporated  
2400 Bernville Road  
Reading, PA 19605 USA

Teleflex®  
www.teleflex.com

LBL010460 (08/20)



REF CV-45854 8.5 Fr. 4 Lumen 20 cm catheter length .032 inch dia. spring-wire guide

# Arrowg+ard Blue Plus® Four-Lumen CVC

