

ARROW



REF EU-04218-N (IPNXXXXXX)

18 Ga. | **1** Lumen | **20** cm catheter length | **.025** inch dia. spring-wire guide

Central Venous Catheterization Set

Contents:

- 1: Indwelling Catheter: 18 Ga. x 20 cm Radiopaque Polyurethane with Integral Suture Wing
- 1: Spring-Wire Guide, Nitinol, Marked: .025" (0.64 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer
- 1: Introducer Needle: 18 Ga. x 2-3/4" (7 cm) XTW
- 1: Syringe: 5 mL Luer-Slip
- 1: Tissue Dilator: 7 Fr. (2.3 mm) x 10.2 cm
- 1: SecondSite™ Adjustable Hub: Catheter Clamp
- 1: SecondSite™ Adjustable Hub: Fastener
- 1: SharpsAway II™ Locking Disposal Cup
- 1: Safety Scalpel: #11


All components are CE 0086 unless otherwise noted.

Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.

Not made with natural rubber latex.

Fluid path components are non-pyrogenic.

- be** Набор для центральной венозной катетеризации
- hr** Centralni venski kateter Komplet
- cs** Centrální žilní katetrizační set
- da** Centralvenøs kateteriserings-sæt
- nl** Centraal veneuze katheterisatieset
- et** Tsentraalveeni kateteriseerimiskomplekt
- fi** Keskuslaskimokatetrisarja
- fr** Jeu de cathétérisme veineux central
- de** Besteck für zentralen Venenkatheter
- el** Σετ Καθετηριασμού Κεντρικής Φλέβας
- hu** Centrális vénás katéterezési szett
- is** Holæðaleggsett
- it** Set per cateterismo della vena centrale
- lt** Centrinės venos kateterizavimo rinkinys
- no** Sentralt venekateteriseringssett
- pl** Zestaw do cewnikowania żył centralnych
- pt** Conjunto de Cateterização Venosa Central
- ro** Set de cateterizare venoasă centrală
- ru** Комплект для установки центрального венозного катетера
- sr** Komplet za Centralnu vensku kateterizaciju
- sk** Centrálna venózna katetrizačná súprava
- sl** Centralni venski komplet za katetrizacijo
- es** Conjunto de cateterización venosa central
- sv** Set för central venkateterisering
- tr** Santral Venöz Kateterizasyon Seti
- uk** Комплект для встановлення центрального венозного катетера



Lumen	Priming Volume* (ml)	Flow Rate† (ml/hr)
Distal (18 Ga.)	0.20	2300

* Priming volumes are approximate and are done without the injection cap. Injection cap priming volume is 0.17 ml.

† Flow rates are done with normal saline, room temperature, 100 cm head height and represent approximate flow capabilities.

Lot Number
 YYYY-MM-DD

(01)XXXXXXXXXXXXXX
 (17)YYMMDD
 (10)LotNumber

Teleflex Medical
 IDA Business and
 Technology Park
 Dublin Road, Athlone
 Co. Westmeath, Ireland

Arrow International, Inc.
 Subsidiary of Teleflex Incorporated
 2400 Bernville Road
 Reading, PA 19605 USA

www.teleflex.com

LBL010479 (02/18)



REF EU-04218-N **18** Ga. | **1** Lumen | **20** cm catheter length | **.025** inch dia. spring-wire guide

Central Venous Catheterization Set

ARROW