

ARROW



REF

EU-25703-EN (IPNXXXXXX)

7 Fr. | **3** Lumen | **20** cm catheter length | **.032** inch dia. spring-wire guide

Multi-Lumen Central Venous Catheterization Set with Blue FlexTip® ARROWg+ard Blue® Catheter

Contents:

- 1: Multi-Lumen Indwelling Catheter: 7 Fr. x 20 cm Radiopaque Polyurethane with Blue FlexTip®, ARROWg+ard® Antimicrobial Surface Treatment[†], Extension Line Clamps
- 1: Spring-Wire Guide, Nitinol, Marked: .032" (0.81 mm) dia. x 23-5/8" (60 cm) (Straight Soft Tip on One End – "J" Tip on Other) with Arrow Advancer with ECG Mark
- 1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW
- 1: Pressure Transduction Probe
- 1: Arrow® Raulerson Spring-Wire Introduction Syringe: 5 mL
- 1: Tissue Dilator: 8.5 Fr. (2.8 mm) x 10.2 cm
- 1: SecondSite™ Adjustable Hub: Fastener
- 1: SecondSite™ Adjustable Hub: Catheter Clamp

[†]Licensed under US Patent No. 6,706,024.


All components are CE 0086 unless otherwise noted.

Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.

Not made with natural rubber latex. Store below 25°C (77°F). Avoid excessive heat above 40°C (104°F).

Fluid path components are non-pyrogenic.

Contraindications: The ARROWg+ard Blue® antimicrobial catheter is contraindicated for patients with known hypersensitivity to chlorhexidine, silver sulfadiazine and/or sulfa drugs (refer to product instructions for references).



Lumen	Priming Volume* (ml)	Flow Rate† (ml/hr)
Distal (16 Ga.)	0.44	3100
Medial (18 Ga.)	0.39	1500
Proximal (18 Ga.)	0.39	1600

* Priming volumes are approximate and are done without the injection cap.

† Flow rates are done with normal saline, room temperature, 100 cm head height and represent approximate flow capabilities.

- be** Многопросветный центральный венозный катетер Набор: Катетер Blue с наконечником Blue FlexTip ARROWg+ard Blue Catheter
- hr** Komplet za višelumensku centralno vensku kateterizaciju sa Blue FlexTip ARROWg+ard Blue kateterom
- cs** Višelumenový centrální žilní katetrizační set s katetrem Blue FlexTip ARROWg+ard Blue
- da** Multilumen centralvenøs kateteriseringssett med Blue FlexTip ARROWg+ard Blue kateter
- nl** Multi-lumen centraal veneuze katheterisatieset met Blue FlexTip ARROWg+ard Blue katheter
- et** Mitme valendikuga tsentraalse veenikateetri komplekt koos Blue FlexTip ARROWg+ard Blue katetriga
- fi** Monilumeninen keskuslaskimokatetrisarja, jossa Blue FlexTip ARROWg+ard Blue-katetri
- fr** Jeu de cathétérisme veineux central à lumières multiples avec cathéter Blue FlexTip ARROWg+ard Blue
- de** Besteck für viellumigen, zentralen Venenkatheter mit Blue-FlexTip ARROWg+ard Blue-Katheter
- el** Σετ Καθητηριασμού Κεντρικής Φλέβας με Πολλαπλούς Αυλούς και Καθετήρα ARROWg+ard Blue με Blue FlexTip
- hu** Többlumenű, centrális vénás katéterezési szett, Blue FlexTip ARROWg+ard Blue katéterrel
- is** Fjöl lúmena holæðaleggsett með Blue FlexTip ARROWg+ard Blue æðalegg
- it** Set per cateterismo della vena centrale multilume con catetere Blue FlexTip ARROWg+ard Blue
- lt** Daugiakanaliscentrinės venos kateterizavimo mini rinkinys su „Blue FlexTip ARROWg+ard Blue“ antimikrobinu kateteriu
- no** Multilumen Sentral venekateteriseringssett med Blue FlexTip ARROWg+ard Blue-kateter
- pl** Zestaw do cewnikowania żyły centralnej z cewnikami wielokanalowymi z końcówką Blue FlexTip i powłoką antybakteryjną ARROWg+ard Blue
- pt** Conjunto de Cateterização Venosa Central com Lúmen Múltiplo e Cateter Blue FlexTip, ARROWg+ard Blue
- ro** Set de cateterizare venoasă centrală multi-lumen cu cateter Blue FlexTip ARROWg+ard Blue
- ru** Набор с многопросветным центральным венозным катетером ARROWg+ard Blue с кончиком Blue FlexTip
- sr** Višelumenski centralni venski kateterizacioni Set sa Blue FlexTip ARROWg+ard Blue kateterom
- sk** Viačlúmenová centrálna venózna katetrizačná súprava s katétrom Blue FlexTip ARROWg+ard Blue
- sl** Večlumenski centralni venski komplet za katetrizacijo s katetrom Blue FlexTip ARROWg+ard Blue
- es** Conjunto de cateterización venosa central de luz múltiple con catéter Blue FlexTip ARROWg+ard Blue
- sv** Set för flerkanalig central venkatetrisering med Blue FlexTip ARROWg+ard Blue-kateter
- tr** ARROWg+ard Blue Multi Lümen Santral Venöz Kateterizasyon Seti, Blue FlexTip Kateter ile
- uk** Багатопрсвітний центральний венозний катетер Набор: катетер Blue с кінчиком Blue FlexTip ARROWg+ard Blue Catheter

LotNumber
 YYYY-MM-DD

(01)XXXXXXXXXXXXXX
 (17)YYMMDD
 (10)LotNumber

Arrow International, Inc.
 Subsidiary of Teleflex Incorporated
 2400 Bernville Road
 Reading, PA 19605 USA

Teleflex
 www.teleflex.com
 LBL010518 (02/18)



REF **EU-25703-EN**

7 Fr. | **3** Lumen | **20** cm catheter length | **.032** inch dia. spring-wire guide

Multi-Lumen Central Venous Catheterization Set with Blue FlexTip® ARROWg+ard Blue® Catheter

ARROW