



REF CS-12402-J (IPNXXXXXX)

4^{Fr.} 2^{Lumen} 5^{cm} catheter length .018^{inch dia. spring-wire guide}

Pediatric Two-Lumen Central Venous Catheterization Set with Blue FlexTip® Catheter


Contents:

- 1: Two-Lumen Indwelling Catheter: 4 Fr. x 2" (5 cm) Radiopaque Polyurethane with Blue FlexTip®, Extension Line Clamps
- 1: Spring-Wire Guide, Marked: .018" (0.46 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer
- 1: Catheter: 22 Ga. x 1" (2.5 cm) Radiopaque over 25 Ga. TW Introducer Needle
- 1: Introducer Needle: 21 Ga. x 1-1/2" (3.81 cm) TW
- 2: Syringe: 3 mL Luer-Slip
- 1: Tissue Dilator: 5 Fr. (1.8 mm) x 7.5 cm
- 1: Tissue Dilator: 5 Fr. (1.8 mm) x 4.1 cm
- 2: Dust Cap: Non-Vented
- 1: SecondSite™ Adjustable Hub: Catheter Clamp
- 1: SecondSite™ Adjustable Hub: Fastener
- 1: Drape: 18" x 26" (46 cm x 66 cm) with 3" (7.6 cm) fenestration
- 1: Scalpel: #11
- 3: Gauze Pad: 2" x 2" (5 cm x 5 cm)
- 5: Gauze Pad: 4" x 4" (10 cm x 10 cm)

Rx only. Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.

Not made with natural rubber latex.

Fluid path components are non-pyrogenic.



Lumen	Priming Volume* (ml)	Flow Rate† (ml/hr)
Distal (22 Ga.)	0.24	930
Proximal (22 Ga.)	0.29	920

* Priming volumes are approximate and are done without accessories.

† Based on 100 cm head height.



LOT LotNumber
YYYY-MM-DD



(01)XXXXXXXXXXXXXX
(17)YYMMDD
(10)LotNumber

Arrow International, Inc.
Subsidiary of Teleflex Incorporated
2400 Bernville Road
Reading, PA 19605 USA



LBL010739 (02/19)



REF CS-12402-J

4^{Fr.} 2^{Lumen} 5^{cm} catheter length .018^{inch dia. spring-wire guide}

Pediatric Two-Lumen Central Venous Catheterization Set with Blue FlexTip® Catheter

