



**REF SA-12853** (IPNXXXXXX)

**8.5** Fr. **3** Lumen **16** cm catheter length **.032** inch dia. spring-wire guide

# Multi-Lumen Central Venous Catheterization Set with Blue FlexTip® Catheter

**Contents:**

- 1: Multi-Lumen Indwelling Catheter: 8.5 Fr. x 16 cm Radiopaque Polyurethane with Blue FlexTip®, Extension Line Clamps
- 1: Spring-Wire Guide, Nitinol, Marked: .032" (0.81 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer
- 1: Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle
- 1: Introducer Needle: Echogenic 18 Ga. x 2-1/2" (6.35 cm) XTW
- 1: Pressure Transduction Probe
- 1: Arrow® Raulerson Spring-Wire Introduction Syringe: 5 mL
- 1: Tissue Dilator: 9 Fr. (3.0 mm) x 10.2 cm
- 1: SecondSite™ Adjustable Hub: Fastener
- 1: SecondSite™ Adjustable Hub: Catheter Clamp
- 1: Drape: 29-1/2" x 35-1/2" (75 cm x 90 cm) with 4" (10 cm) fenestration, with adhesive
- 2: Drape: 19-5/8" x 29-1/2" (50 cm x 75 cm) with adhesive
- 1: Safety Scalpel: #11
- 2: Prep Sponge Swab
- 3: Valve: MicroCLAVE®<sup>1</sup> Neutral Displacement Connector


<sup>1</sup>A registered trademark of ICU Medical, Inc.

All components are CE 0086 unless otherwise noted.

**Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.**

Not made with natural rubber latex.

Fluid path components are non-pyrogenic.



**LOT** LotNumber  
YYYY-MM-DD



(01)XXXXXXXXXXXXXXXX  
(17)YYMMDD  
(10)LotNumber



Arrow International, Inc.  
Subsidiary of Teleflex Incorporated  
2400 Bernville Road  
Reading, PA 19605 USA

**Teleflex®**  
www.teleflex.com  
LBL010767 (08/19)



**REF SA-12853**

**8.5** Fr. **3** Lumen **16** cm catheter length **.032** inch dia. spring-wire guide

# Multi-Lumen Central Venous Catheterization Set with Blue FlexTip® Catheter

