



REF SA-22855-EU (IPNXXXXXX)

ARROW®

8.5 Fr. 5 Lumen 16 cm catheter length .032 inch dia. spring-wire guide

Multi-Lumen Central Venous Catheterization Kit with Blue FlexTip® ARROWg+ard Blue® Catheter

Contents:

- 1: Multi-Lumen Indwelling Catheter: 8.5 Fr. x 16 cm Radiopaque Polyurethane with Blue FlexTip®, ARROWg+ard® Antimicrobial Surface Treatment¹, Extension Line Clamps
- 1: Spring-Wire Guide, Nitinol, Marked: .032" (0.81 mm) dia. x 23-5/8" (60 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer with ECG Mark
- 1: Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle
- 1: Introducer Needle: Echogenic 18 Ga. x 2-1/2" (6.35 cm) XTW
- 1: Pressure Transduction Probe
- 1: Arrow® Raulerson Spring-Wire Introduction Syringe: 5 mL
- 1: Tissue Dilator: 9 Fr. (3.0 mm) x 10.2 cm
- 1: SecondSite™ Adjustable Hub: Fastener
- 1: SecondSite™ Adjustable Hub: Catheter Clamp
- 1: Drape: 29-1/2" x 35-1/2" (75 cm x 90 cm) with 4" (10 cm) fenestration, with adhesive
- 2: Drape: 19-5/8" x 29-1/2" (50 cm x 75 cm) with adhesive
- 1: Safety Scalpel: #11
- 2: Prep Sponge Swab
- 5: Valve: MicroCLAVE®² Neutral Displacement Connector


¹Licensed under US Patent No. 6,706,024.
²A registered trademark of ICU Medical, Inc.

All components are CE 0086 unless otherwise noted.

Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.

Not made with natural rubber latex.
Store below 25°C (77°F). Avoid excessive heat above 40°C (104°F).
Fluid path components are non-pyrogenic.

Contraindications: The Arrowg+ard Blue antimicrobial catheter is contraindicated for patients with known hypersensitivity to chlorhexidine, silver sulfadiazine and/or sulfa drugs (refer to product instructions for references).



Lumen	Priming Volume* (mL)	Flow Rate† (mL/hr)
Distal (16 Ga.)	0.58	3650
Medial 1 (14 Ga.)	0.71	6590
Medial 2 (18 Ga.)	0.36	1260
Medial 3 (18 Ga.)	0.38	1430
Proximal (18 Ga.)	0.4	1460

* Priming volumes are approximate and are done without accessories.
† Based on 100 cm head height.



LOT LotNumber
YYYY-MM-DD



(01)XXXXXXXXXXXXXXXX
(17)YYMMDD
(10)LotNumber



Arrow International, Inc.
Subsidiary of Teleflex Incorporated
2400 Bernville Road
Reading, PA 19605 USA

Teleflex®
www.teleflex.com
LBL010809 (08/19)



REF SA-22855-EU

8.5 Fr. 5 Lumen 16 cm catheter length .032 inch dia. spring-wire guide

Multi-Lumen Central Venous Catheterization Kit with Blue FlexTip® ARROWg+ard Blue® Catheter

ARROW®