

Dimensions:
Width: 6-9/64" (15.6 cm)
Length: 8-21/32" (22 cm)

PB-50014-100B
PBZ-50014-103A

T50014100A6 Dimensions: 5 7/8" x 7 5/8"

ARROW®



REF CV-12703-E (IPNXXXXXX)

7 Fr. 3 Lumen 16 cm catheter length .032 inch dia. spring-wire guide

Multi-Lumen Central Venous Catheterization Set with Blue FlexTip® Catheter

Contents:

- 1: Multi-Lumen Indwelling Catheter: 7 Fr. x 16 cm Radiopaque Polyurethane with Blue FlexTip®, Extension Line Clamps
- 1: Spring-Wire Guide, Marked: .032" (0.81 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer
- 1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW
- 1: Pressure Transduction Probe
- 1: Arrow® Raulerson Spring-Wire Introduction Syringe: 5 mL
- 1: Tissue Dilator: 8.5 Fr. (2.8 mm) x 10.2 cm
- 3: Dust Cap: Non-Vented
- 1: SecondSite™ Adjustable Hub: Fastener
- 1: SecondSite™ Adjustable Hub: Catheter Clamp

All components are CE 0086 unless otherwise noted.

Rx only. **Warning:** Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.

California Prop. 65


⚠ **WARNING:** Cancer and Reproductive Harm.

www.P65Warnings.ca.gov

Not made with natural rubber latex.

Fluid path components are non-pyrogenic.

- fr** Jeu de cathétérisme veineux central à lumières multiples avec cathéter Blue FlexTip
- de** Besteck für viellumigen, zentralen Venenkatheter mit Blue FlexTip Katheter
- it** Set per cateterismo della vena centrale multilume con catetere Blue FlexTip
- pl** Zestaw do cewnikowania żyły centralnej z cewnikami wielokanałowymi z końcówką Blue FlexTip
- pt** Conjunto de Cateterização Venosa Central com Lúmen Múltiplo e Cateter Blue FlexTip
- ru** Набор с многопросветным центральным венозным катетером с кончиком Blue FlexTip
- sl** Večlumenski centralni venski komplet za katetizacijo katetrom Blue FlexTip
- es** Conjunto de cateterización venosa central de luz múltiple con catéter Blue FlexTip
- sv** Set för flerkanalig central venkatetrering med Blue FlexTip kateter
- tr** Multi Lúmen Santral Venöz Kateterizasyon Seti, Mavi esnek uçlu Kateter ile



| Lumen | Priming Volume* (ml) | Flow Rate† (ml/hr) |
|-------------------|----------------------|--------------------|
| Distal (16 Ga.) | 0.39 | 3400 |
| Medial (18 Ga.) | 0.35 | 1800 |
| Proximal (18 Ga.) | 0.37 | 1900 |

* Priming volumes are approximate and are done without accessories.

† Based on 100 cm head height.

Catheter Cross Section



LOT LotNumber
YYYY-MM-DD



(01)XXXXXXXXXXXXXX
(17)YYMMDD
(10)LotNumber

CE 0086 Teleflex Medical
IDA Business and
Technology Park
Dublin Road, Athlone
Co. Westmeath, Ireland

Arrow International, Inc.
Subsidiary of Teleflex Incorporated
2400 Bernville Road
Reading, PA 19605 USA

teleflex®
www.teleflex.com

LBL010824 (01/19)



7 Fr. 3 Lumen 16 cm catheter length .032 inch dia. spring-wire guide

REF CV-12703-E

Multi-Lumen Central Venous Catheterization Set with Blue FlexTip® Catheter

ARROW®