

Dimensions:
Width: 6-9/64" (15.6 cm)
Length: 8-21/32" (22 cm)

PB-50014-100B
PBZ-50014-103A

T50014100A6 Dimensions: 5 7/8" x 7 5/8"

ARROW®



REF CV-16702-E (IPNXXXXXX)

7 Fr. 2 Lumen 20 cm catheter length .032 inch dia. spring-wire guide

Two-Lumen Central Venous Catheterization Set with Blue FlexTip® Catheter

Contents:

- 1: Two-Lumen Indwelling Catheter: 7 Fr. x 20 cm Radiopaque Polyurethane with Blue FlexTip®, Extension Line Clamps
- 1: Spring-Wire Guide, Marked: .032" (0.81 mm) dia. x 23-5/8" (60 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer
- 1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW
- 1: Pressure Transduction Probe
- 1: Arrow® Raulerson Spring-Wire Introduction Syringe: 5 mL
- 1: Tissue Dilator: 8.5 Fr. (2.8 mm) x 10.2 cm
- 2: Dust Cap: Non-Vented
- 1: SecondSite™ Adjustable Hub: Fastener
- 1: SecondSite™ Adjustable Hub: Catheter Clamp

All components are CE 0086 unless otherwise noted.

Rx only. **Warning:** Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.

California Prop. 65


⚠ **WARNING:** Cancer and Reproductive Harm.

www.P65Warnings.ca.gov

Not made with natural rubber latex.

Fluid path components are non-pyrogenic.

- fr** Jeu de cathétérisme veineux central à deux lumières avec cathéter Blue FlexTip
- de** Besteck für zweilumigen, zentralen Venenkatheter mit Blue FlexTip Katheter
- it** Set per cateterismo della vena centrale a due lumi con catetere Blue FlexTip
- pl** Zestaw do cewnikowania żyły centralnej z cewnikiem 2-kanalowym z końcówką Blue FlexTip
- pt** Conjunto de Cateterização Venosa Central com Lúmen Duplo e Cateter Blue FlexTip
- ru** Двухпросветный центральный венозный катетер с наконечником Blue FlexTip
- sl** Dvolumenski centralni venski pripor za katetrizacijo s katetrom Blue FlexTip
- es** Conjunto de cateterización venosa central de dos lúmenes con catéter Blue FlexTip
- sv** Set för tvåkanalig central venkatetrisering med Blue FlexTip kateter
- tr** İki Lümenli Santral Venöz Kateterizasyon Seti, Mavi esnek uçlu kateterli



Lumen	Priming Volume* (ml)	Flow Rate† (ml/hr)
Distal (16 Ga.)	0.47	3700
Proximal (16 Ga.)	0.45	3000

* Priming volumes are approximate and are done without accessories.
† Based on 100 cm head height.

Catheter Cross Section



LOT LotNumber
YYYY-MM-DD



(01)XXXXXXXXXXXXXX
(17)YYMMDD
(10)LotNumber

CE REP Teleflex Medical
IDA Business and
Technology Park
0086 Dublin Road, Athlone
Co. Westmeath, Ireland

Arrow International, Inc.
Subsidiary of Teleflex Incorporated
2400 Bernville Road
Reading, PA 19605 USA

teleflex®
www.teleflex.com

LBL010828 (01/19)



7 Fr. 2 Lumen 20 cm catheter length .032 inch dia. spring-wire guide

REF CV-16702-E

Two-Lumen Central Venous Catheterization Set with Blue FlexTip® Catheter

ARROW®