

Dimensions:  
Width: 6-9/64" (15,6 cm)  
Length: 8-21/32" (22 cm)

PB-50014-100B  
PBZ-50014-103A

T50014100A6 Dimensions: 5 7/8" x 7 5/8"

ARROW®



REF CV-24703-E (IPNXXXXXX)

7 Fr. 3 Lumen 30 cm catheter length .032 inch dia. spring-wire guide

## Multi-Lumen Central Venous Catheterization Set with Blue FlexTip® ARROWg+ard Blue® Catheter

### Contents:

- 1: Multi-Lumen Indwelling Catheter: 7 Fr. x 30 cm Radiopaque Polyurethane with Blue FlexTip®, ARROWg+ard® Antimicrobial Surface Treatment<sup>1</sup>, Extension Line Clamps
- 1: Spring-Wire Guide, Marked: .032" (0.81 mm) dia. x 26-3/4" (68 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer with ECG Mark
- 1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW
- 1: Pressure Transduction Probe
- 1: Arrow® Raulerson Spring-Wire Introduction Syringe: 5 mL
- 1: Tissue Dilator: 8.5 Fr. (2.8 mm) x 10.2 cm
- 3: Dust Cap: Non-Vented
- 1: SecondSite™ Adjustable Hub: Fastener
- 1: SecondSite™ Adjustable Hub: Catheter Clamp

<sup>1</sup>Licensed under US Patent No. 6,706,024.

All components are CE 0086 unless otherwise noted.

- fr** Jeu de cathétérisme veineux central à lumières multiples avec cathéter Blue FlexTip ARROWg+ard Blue
- de** Besteck für viellumigen, zentralen Venenkatheter mit Blue FlexTip ARROWg+ard Blue Katheter
- it** Set per cateterismo della vena centrale multilume con catetere Blue FlexTip ARROWg+ard Blue
- pl** Zestaw do cewnikowania żyły centralnej z cewnikiem wielokanałowym z końcówką Blue FlexTip i powłoką antybakteryjną ARROWg+ard Blue
- pt** Conjunto de Cateterização Venosa Central com Lúmen Múltiplo e Cateter Blue FlexTip, ARROWg+ard Blue
- ru** Набор с многопросветным центральным венозным катетером ARROWg+ard Blue с кончиком Blue FlexTip
- sl** Večlumenski centralni venški komplet za katetizacijo s katetrom Blue FlexTip ARROWg+ard Blue
- es** Conjunto de cateterización venosa central de luz múltiple con catéter Blue FlexTip ARROWg+ard Blue
- sv** Set för flerkanalig central venkatetrering med Blue FlexTip ARROWg+ard Blue kateter
- tr** Multi Lúmen Santral Venöz Kateterizasyon Seti, Mavi esnek uçlu ARROWg+ard Blue Kateterli

Rx only. **Warning:** Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.


California Prop. 65

**WARNING:** Cancer and Reproductive Harm.  
www.P65Warnings.ca.gov

Not made with natural rubber latex. Store below 25°C (77°F). Avoid excessive heat above 40°C (104°F).

Fluid path components are non-pyrogenic.

**Contraindications:** The Arrowg+ard Blue antimicrobial catheter is contraindicated for patients with known hypersensitivity to chlorhexidine, silver sulfadiazine and/or sulfa drugs (refer to product instructions for references).



Lumen	Priming Volume* (ml)	Flow Rate† (ml/hr)
Distal (16 Ga.)	0.55	2400
Medial (18 Ga.)	0.45	1070
Proximal (18 Ga.)	0.5	1160

\* Priming volumes are approximate and are done without accessories.

† Based on 100 cm head height.

Catheter Cross Section



**LOT** LotNumber  
YYYY-MM-DD



(01)XXXXXXXXXXXXXX  
(17)YYMMDD  
(10)LotNumber

**CE** Teleflex Medical  
IDA Business and  
Technology Park  
0086 Dublin Road, Athlone  
Co. Westmeath, Ireland

Arrow International, Inc.  
Subsidiary of Teleflex Incorporated  
2400 Bernville Road  
Reading, PA 19605 USA

**teleflex®**  
www.teleflex.com

LBL010834 (02/19)



7 Fr. 3 Lumen 30 cm catheter length .032 inch dia. spring-wire guide

REF CV-24703-E

Multi-Lumen Central Venous Catheterization Set with Blue FlexTip® ARROWg+ard Blue® Catheter

ARROW®