

Dimensions:
Width: 10"
Length: 17-3/4"

PB-05000-111A

T-05600-102S-04 Dimensions: 9 7/8" x 15 7/8"



REF ASK-04701-KR (IPNXXXXXX)

14 Ga. **1** Lumen **20** cm catheter length **.032** inch dia. spring-wire guide

Single-Lumen CVC Kit with Blue FlexTip® Catheter

Contents:

- 1: Indwelling Catheter: 14 Ga. x 20 cm Radiopaque Polyurethane with Blue FlexTip®, Integral Extension Line, Extension Line Clamp, Side Holes, Integral Suture Wing, Injection Site Cap
- 1: Spring-Wire Guide, Marked: .032" (0.81 mm) dia. x 23-5/8" (60 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer
- 1: Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle
- 1: Injection Needle: 22 Ga. x 1-1/2" (3.81 cm)
- 1: Injection Needle: 25 Ga. x 1-1/2" (3.81 cm)
- 1: Injection Needle: 18 Ga. x 1-1/2" (3.81 cm)
- 1: Introducer Needle: 18 Ga. x 1-1/2" (3.81 cm) XTW
- 1: Introducer Needle: Echogenic 18 Ga. x 2-1/2" (6.35 cm) XTW
- 1: Pressure Transduction Probe
- 1: Introducer Needle: 22 Ga. x 2-3/4" (7 cm) TW
- 1: Syringe: 10 mL Luer-Slip
- 1: Syringe: 5 mL Luer-Slip
- 1: Arrow® Raulerson Spring-Wire Introduction Syringe: 5 mL
- 1: Tissue Dilator: 7 Fr. (2.3 mm) x 10.2 cm
- 1: SecondSite™ Adjustable Hub: Fastener
- 1: SecondSite™ Adjustable Hub: Catheter Clamp
- 2: Medicine Cup: 60 mL
- 1: SharpsAway® Disposal Cup
- 1: Drape: 35" x 50" (89 cm x 127 cm) with 4" (10 cm) fenestration, with adhesive
- 1: Drape: 30" x 44" (76 cm x 112 cm)
- 1: Needle Holder
- 1: Scalpel: #11
- 1: Stopcock: 4-way
- 2: Prep Sponge Swab
- 6: Gauze Pad: 4" x 4" (10 cm x 10 cm)
- 1: Dressing: OpSite®¹ 4-3/4" x 4" (12 cm x 10 cm)
- 1: Suture: 3-0 Silk with Curved Needle
- 2: Tray: Prep

¹A registered trademark of T.J. Smith & Nephew, Limited.

FOR EXPORT ONLY


Shonin Approval Number 20200BZY00019000

Rx only.

Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.

Not made with natural rubber latex.

Fluid path components are non-pyrogenic.



LOT LotNumber
YYYY-MM-DD



(01)XXXXXXXXXXXXXX
(17)YYMMDD
(10)LotNumber

Arrow International, Inc.
Subsidiary of Teleflex Incorporated
2400 Bernville Road
Reading, PA 19605 USA

Teleflex®
www.teleflex.com
LBL010999 (12/18)



REF ASK-04701-KR

14 Ga. **1** Lumen **20** cm catheter length **.032** inch dia. spring-wire guide

Single-Lumen CVC Kit with Blue FlexTip® Catheter

