

Dimensions:
Width: 10"
Length: 17-3/4"

PB-05000-111A

T-05600-102S-04 Dimensions: 9 7/8" x 15 7/8"



REF ASK-14854-KR (IPNXXXXXX)

8.5 Fr. **4** Lumen **30** cm catheter length **.032** inch dia. spring-wire guide

Quad-Lumen CVC Kit with Blue FlexTip® Catheter

Contents:

- 1: Quad-Lumen Indwelling Catheter: 8.5 Fr. x 30 cm Radiopaque Polyurethane with Blue FlexTip®, Extension Line Clamps
- 1: Spring-Wire Guide, Marked: .032" (0.81 mm) dia. x 26-3/4" (68 cm) (Straight Soft Tip on One End – "J" Tip on Other) with Arrow Advancer
- 1: Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle
- 1: Injection Needle: 22 Ga. x 1-1/2" (3.81 cm)
- 1: Injection Needle: 25 Ga. x 1-1/2" (3.81 cm)
- 1: Injection Needle: 18 Ga. x 1-1/2" (3.81 cm)
- 1: Introducer Needle: 18 Ga. x 1-1/2" (3.81 cm) XTW
- 1: Introducer Needle: Echogenic 18 Ga. x 2-1/2" (6.35 cm) XTW
- 1: Pressure Transduction Probe
- 1: Introducer Needle: 22 Ga. x 2-3/4" (7 cm) TW
- 1: Syringe: 10 mL Luer-Slip
- 1: Syringe: 5 mL Luer-Slip
- 1: Arrow® Raulerson Spring-Wire Introduction Syringe: 5 mL
- 1: Tissue Dilator: 8.5 Fr. (2.8 mm) x 10.2 cm
- 4: Dust Cap: Non-Vented
- 1: SecondSite™ Adjustable Hub: Fastener
- 1: SecondSite™ Adjustable Hub: Catheter Clamp
- 2: Medicine Cup: 60 mL
- 1: SharpsAway® Disposal Cup
- 1: Drape: 35" x 50" (89 cm x 127 cm) with 4" (10 cm) fenestration, with adhesive
- 1: Drape: 30" x 44" (76 cm x 112 cm)
- 1: Needle Holder
- 1: Scalpel: #11
- 4: Stopcock: 4-way
- 2: Prep Sponge Swab
- 6: Gauze Pad: 4" x 4" (10 cm x 10 cm)
- 1: Dressing: OpSite® 4-3/4" x 4" (12 cm x 10 cm)
- 1: Suture: 3-0 Silk with Curved Needle
- 2: Tray: Prep


*A registered trademark of T.J. Smith & Nephew, Limited.
FOR EXPORT ONLY
Shonin Approval Number 20200BZY00019000

Rx only.

Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.

Not made with natural rubber latex.

Fluid path components are non-pyrogenic.



Lumen	Priming Volume* (mL)	Flow Rate† (mL/hr)
Distal (16 Ga.)	0.52	2100
Medial 1 (14 Ga.)	0.66	3900
Medial 2 (18 Ga.)	0.44	1000
Proximal (18 Ga.)	0.47	1000

* Priming volumes are approximate and are done without accessories.
† Based on 100 cm head height.



LOT LotNumber
YYYY-MM-DD



(01)XXXXXXXXXXXXXX
(17)YYMMDD
(10)LotNumber

Arrow International, Inc.
Subsidiary of Teleflex Incorporated
2400 Bernville Road
Reading, PA 19605 USA

Teleflex®
www.teleflex.com
LBL011008 (03/20)



REF ASK-14854-KR

8.5 Fr. **4** Lumen **30** cm catheter length **.032** inch dia. spring-wire guide

Quad-Lumen CVC Kit with Blue FlexTip® Catheter

