

Dimensions:  
Width: 8-3/4" (22,2 cm)  
Length: 11-1/8" (28,3 cm)

PB-10001-000A  
PBZ-10001-000A

T12703100B2 Dimensions:  
Width: 8.5"  
Length: 10.125

ARROW®



REF EU-15802-CVTS (IPNXXXXXX)

8 Fr. | 2 Lumen | 20 cm catheter length | .032 inch dia. spring-wire guide

## Two-Lumen Central Venous Catheterization Set with Blue FlexTip® Catheter

### Contents:

- 1: Two-Lumen Indwelling Catheter: 8 Fr. x 20 cm Radiopaque Polyurethane with Blue FlexTip®, Extension Line Clamps
- 1: Spring-Wire Guide, Marked: .032" (0.81 mm) dia. x 23-5/8" (60 cm) (Straight Soft Tip on One End – "J" Tip on Other) with Arrow Advancer with ECG Mark
- 1: Introducer Needle: Echogenic 18 Ga. x 2-1/2" (6.35 cm) XTW
- 1: Syringe: 5 mL Luer-Slip
- 1: Tissue Dilator: 9 Fr. (3.0 mm) x 10.2 cm
- 1: ECG Extension Cable
- 2: Injection Site Cap
- 1: SecondSite™ Adjustable Hub: Fastener
- 1: SecondSite™ Adjustable Hub: Catheter Clamp
- 1: SharpsAway II™ Locking Disposal Cup
- 1: Safety Scalpel: #11


All components are CE 0086 unless otherwise noted.

**Warning:** Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.

Not made with natural rubber latex.

Fluid path components are non-pyrogenic.

- be** Двухпросветный центральный венозный катетер с наконечником Blue FlexTip
- hr** Dvolumenski komplet za centralnu vensku kateterizaciju s plavim savitljivim vrhom FlexTip
- cs** Dvolumenový centrální žilní katetrizační set s katetrem Blue FlexTip
- da** To-lumen centralvenekateteriseringsset med Blue FlexTip kateter
- nl** Centraal veneuze twee-lumen katheterisatieset met Blue FlexTip katheter
- et** Kahe luumeniga tsentraalveeni kateteriseerimiskomplekt sinise painduva otsaga Blue kateetriiga
- fi** Kaksiluumen-keskuslaskimokatetri setti, jossa Blue FlexTip katetri
- fr** Jeu de cathétérisme veineux central à deux lumières avec cathéter Blue FlexTip
- de** Besteck für zweilumigen, zentralen Venenkatheter mit Blue FlexTip Katheter
- el** Σετ Καθετηριασμού Κεντρικής Φλέβας Διπλού Αυλού με Καθετήρα Blue FlexTip
- hu** Két lumenű centrális vénás katéterezési készlet Blue FlexTip katéterrel
- is** Tveggja-lúmena sett fyrir holæðalegg með Blue FlexTip æðalegg
- it** Set per cateterismo della vena centrale a due lumi con catetere Blue FlexTip
- lt** Dviejų kanalų centrinės venos kateterizavimo rinkinys su „Blue FlexTip“ kateteriu
- no** To-Lumen Sentraalt venekateteriseringssett med Blue FlexTip kateter
- pl** Zestaw do cewnikowania żył centralnych, z dwukanałowym cewnikiem z końcówką Blue FlexTip i powłoką
- pt** Conjunto de Cateterização Venosa Central com Lúmen Duplo e Cateter Blue FlexTip
- ro** Trusa de cateterizare venoasă centrală cu două lumene, cu cateter Blue FlexTip
- ru** Двухпросветный центральный венозный катетер с наконечником Blue FlexTip
- sr** Dvo-lumenski centralni venski kateterizacioni Set sa Blue FlexTip kateterom
- sk** Dvojlúmenová katetrizačná súprava na katetrizáciu centrálnej žily s katétrom Blue FlexTip
- sl** Dvolumenski centralni venski pribor za katetrizacijo s katetrom Blue FlexTip
- es** Conjunto de cateterización venosa central de dos luces con catéter Blue FlexTip
- sv** Set för tvåkanalyg central venkateterisering med Blue FlexTip kateter
- tr** Blue FlexTip Kateterli İki Lümenli Santral Venöz Kateterizasyon Seti
- uk** Двохпросвітний центральний венозний катетер з наконечником Blue FlexTip



Lumen	Priming Volume* (ml)	Flow Rate† (ml/hr)
Distal (14 Ga.)	0.75	5600
Proximal (14 Ga.)	0.80	5100

\* Priming volumes are approximate and are done without the injection cap. Injection cap priming volume is 0.17 ml.

† Flow rates are done with normal saline, room temperature, 100 cm head height and represent approximate flow capabilities.



**LOT** LotNumber  
YYYY-MM-DD



(01)XXXXXXXXXXXXXXXXXX  
(17)YYMMDD  
(10)LotNumber

**CE** Teleflex Medical  
IDA Business and  
Technology Park  
Dublin Road, Athlone  
Co. Westmeath, Ireland  
0086

Arrow International, Inc.  
Subsidiary of Teleflex Incorporated  
2400 Bernville Road  
Reading, PA 19605 USA

**Teleflex®**  
www.teleflex.com  
LBL011110 (04/18)



REF EU-15802-CVTS

8 Fr. | 2 Lumen | 20 cm catheter length | .032 inch dia. spring-wire guide

Two-Lumen Central Venous Catheterization Set with Blue FlexTip® Catheter

ARROW