



REF ES-24402 (IPNXXXXXX)

4 Fr. 2 Lumen 13 cm catheter length .018 inch dia. spring-wire guide

Pediatric Arrowgard Blue® Two-Lumen CVC

Contents:

- 1: Arrowgard Blue® Two-Lumen Catheter: 4 Fr. (1.5 mm OD) x 13 cm
- 1: Spring-Wire Guide: .018" (0.46 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on Both Ends)
- 1: Introducer Needle: 22 Ga. x 1-1/16" (2.67 cm) TW
- 1: Syringe: 3 mL Luer-Slip
- 1: Tissue Dilator: 5 Fr. (1.8 mm) x 7.5 cm
- 1: Tissue Dilator: 5 Fr. (1.8 mm) x 4.1 cm
- 2: Dust Cap: Non-Vented
- 1: SecondSite™ Adjustable Hub: Catheter Clamp
- 1: SecondSite™ Adjustable Hub: Fastener

All components are CE 0086 unless otherwise noted.

Rx only Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.

Not made with natural rubber latex. Store below 25°C (77°F). Avoid excessive heat above 40°C (104°F).

Fluid path components are non-pyrogenic.

Contraindications:

The Arrowgard Blue antimicrobial catheter is contraindicated for patients with known hypersensitivity to chlorhexidine, silver sulfadiazine and/or sulfa drugs.

- fr** CVC double lumière Arrowgard Blue pédiatrique
- de** Pädiatrischer Arrowgard Blue zweilumiger ZVK
- it** CVC a due lumi Arrowgard Blue pediatrico
- pl** Dwukanałowy cewnik pediatryczny CVC Arrowgard Blue
- pt** CVC de dois lúmenes Arrowgard Blue pediátrico
- ru** Педиатрический двухпросветный центральный венозный катетер Arrowgard Blue
- sl** Pediatrični centralni venski kateter z dvema svetlinama Arrowgard Blue
- es** CVC de dos luces Arrowgard Blue pediátrico
- sv** Pediatrisk Arrowgard Blue tvålumen-CVK
- tr** Pediatrik Arrowgard Blue İki Lümenli SVK



Lumen	Priming Volume* (mL)	Flow Rate† (mL/hr)
Distal (22 Ga.)	0.26	640
Proximal (22 Ga.)	0.30	290

* Priming volumes are approximate and are done without accessories.
 † Based on 100 cm head height.

LOT LotNumber
 YYYY-MM-DD

(01)XXXXXXXXXXXXXX
 (17)YYMMDD
 (10)LotNumber

CE 0086

EC REP Teleflex Medical
 IDA Business and
 Technology Park
 Dublin Road, Athlone
 Co. Westmeath, Ireland

Arrow International, Inc.
 Subsidiary of Teleflex Incorporated
 2400 Bernville Road
 Reading, PA 19605 USA

Teleflex®
www.teleflex.com
 LBL011175 (04/21)



REF ES-24402

4 Fr. 2 Lumen 13 cm catheter length .018 inch dia. spring-wire guide

Pediatric Arrowgard Blue® Two-Lumen CVC

