

Dimensions:
Width: 6-9/64" (15.6 cm)
Length: 8-21/32" (22 cm)

PB-50014-100B
PBZ-50014-103A

T50014100A6 Dimensions: 5 7/8" x 7 5/8"

ARROW®



REF CV-22702-E (IPNXXXXXX)

7 Fr. 2 Lumen 16 cm catheter length .025 inch dia. spring-wire guide

Two-Lumen Central Venous Catheterization Set with Blue FlexTip® ARROWg+ard Blue® Catheter

Contents:

- 1: Two-Lumen Indwelling Catheter: 7 Fr. x 16 cm Radiopaque Polyurethane with Blue FlexTip®, ARROWg+ard® Antimicrobial Surface Treatment¹, Extension Line Clamps
- 1: Spring-Wire Guide, Marked: .025" (0.64 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer
- 1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) RW
- 1: Pressure Transduction Probe
- 1: Arrow® Raulerson Spring-Wire Introduction Syringe: 5 mL
- 1: Tissue Dilator: 8.5 Fr. (2.8 mm) x 10.2 cm
- 2: Dust Cap: Non-Vented
- 1: SecondSite™ Adjustable Hub: Fastener
- 1: SecondSite™ Adjustable Hub: Catheter Clamp

¹Licensed under US Patent No. 6,706,024.

All components are CE 0086 unless otherwise noted.

- fr** Jeu de cathétérisme veineux central à deux lumières avec cathéter Blue FlexTip ARROWg+ard Blue
- de** Besteck für zweilumigen, zentralen Venenkatheter mit Blue FlexTip ARROWg+ard Blue Katheter
- it** Set per cateterismo della vena centrale a due lumi con catetere Blue FlexTip ARROWg+ard Blue
- pl** Zestaw do cewnikowania żyły centralnej, z cewnikiem 2-kanalowym z końcówką Blue FlexTip i powłoką ARROWg+ard Blue
- pt** Conjunto de Cateterização Venosa Central com Lúmen Duplo e Cateter Blue FlexTip, ARROWg+ard Blue
- ru** Двухпросветный центральный венозный катетер с наконечником Blue FlexTip ARROWg+ard
- sl** Dvolumenski centralni venski pribor za katetrizacijo s katetrom Blue FlexTip ARROWg+ard Blue
- es** Conjunto de cateterización venosa central de dos luces con catéter Blue FlexTip ARROWg+ard Blue
- sv** Set för tvåkanalig central venkatetrisering med Blue FlexTip ARROWg+ard Blue kateter
- tr** İki Lümenli Santral Venöz Kateterizasyon Seti, Mavi esnek uçlu ARROWg+ard Blue Kateterli

Rx only. Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.

California Prop. 65

⚠ **WARNING:** Cancer and Reproductive Harm.
www.P65Warnings.ca.gov


Not made with natural rubber latex. Store below 25°C (77°F). Avoid excessive heat above 40°C (104°F).

Fluid path components are non-pyrogenic.

Contraindications: The Arrowg+ard Blue antimicrobial catheter is contraindicated for patients with known hypersensitivity to chlorhexidine, silver sulfadiazine and/or sulfa drugs (refer to product instructions for references).

Lumen	Priming Volume* (ml)	Flow Rate† (ml/hr)
	Distal (18 Ga.)	0.32
Proximal (14 Ga.)	0.47	5000

* Priming volumes are approximate and are done without accessories.
† Based on 100 cm head height.



Catheter Cross Section



LOT LotNumber
YYYY-MM-DD



(01)XXXXXXXXXXXXXX
(17)YYMMDD
(10)LotNumber

CE Teleflex Medical
IDA Business and
Technology Park
0086 Dublin Road, Athlone
Co. Westmeath, Ireland

Arrow International, Inc.
Subsidiary of Teleflex Incorporated
2400 Bernville Road
Reading, PA 19605 USA

teleflex®
www.teleflex.com

LBL011269 (01/19)



REF CV-22702-E

7 Fr. 2 Lumen 16 cm catheter length .025 inch dia. spring-wire guide

Two-Lumen Central Venous Catheterization Set with Blue FlexTip® ARROWg+ard Blue® Catheter

ARROW®