



**REF** AD-09903 (IPNXXXXXX)

**ARROW**

**9** Fr. | **10** cm sheath length | **3.9** mm sheath O.D. | **.035** inch dia. spring-wire guide

## Percutaneous Sheath Introducer Set with Integral Hemostasis Valve/Side Port for use with 7 - 8 Fr. Catheters

**Contents:**

- 1: Sheath: 9 Fr. x 4" (10 cm) Arrow-Flex® Radiopaque Polyurethane with Integral Hemostasis Valve/Side Port and Tissue Dilator
- 1: Spring-Wire Guide, Marked: .035" (0.89 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End – "J" Tip on Other) with Arrow Advancer
- 1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW
- 1: Syringe: 5 mL Luer-Slip
- 1: Obturator: 8 Fr.
- 1: Dust Cap: Non-Vented
- 1: Stopcock: 4-way High-Flow

**Rx only. Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.**

Not made with natural rubber latex.

Fluid path components are non-pyrogenic.

- fr** Jeu d'introduction de gaines par voie percutanée avec valve hémostatique/port latéral intégrés pour emploi avec les cathéters 7 - 8 Fr.
- de** Perkutanes Schleuseneinführbesteck mit eingebauter Hämostaseventil/Seitenanschluß-Baugruppe zum Gebrauch mit Kathetern der Größen 7 - 8 Fr.
- it** Set per introduzione di guaina percutanea con valvola emostatica/bocchetta laterale incorporate per l'uso con cateteri da 7 - 8 Fr.
- pt** Conjunto de Introdução de Bainha Percutânea com Dispositivo de Válvula de Hemostase Integral e Porta Lateral, para utilizar com Cateteres 7 - 8 Fr.
- es** Conjunto introductor de vaina percutánea con válvula hemostática y orificio lateral integrales para uso con catéteres de 7 - 8 Fr.
- sv** Set för perkutan skyddsintroducerare med integral hemostasventil/sidoöppning för bruk med katetrar mellan 7 - 8 Fr.



**LOT** LotNumber  
YYYY-MM-DD



(01)XXXXXXXXXXXXXX  
(17)YYMMDD  
(10)LotNumber

**EC REP** Teleflex Medical  
IDA Business and  
Technology Park  
Dublin Road, Athlone  
0086 Co. Westmeath, Ireland

Arrow International, Inc.  
Subsidiary of Teleflex Incorporated  
2400 Bernville Road  
Reading, PA 19605 USA

**Teleflex**  
www.teleflex.com  
LBL011296 (06/18)



**REF** AD-09903

**9** Fr. | **10** cm sheath length | **3.9** mm sheath O.D. | **.035** inch dia. spring-wire guide

## Percutaneous Sheath Introducer Set with Integral Hemostasis Valve/Side Port for use with 7 - 8 Fr. Catheters

**ARROW**