

Special Order Product



REF AD-09903-S (IPNXXXXXX)



9 Fr. | **10** cm sheath length | **3.9** mm sheath O.D.

Percutaneous Sheath Introducer Set with Integral Hemostasis Valve/Side Port for use with 7 - 8 Fr. Catheters

Contents:

- 1: Sheath: 9 Fr. x 4" (10 cm) Arrow-Flex® Radiopaque Polyurethane with Integral Hemostasis Valve/Side Port and Tissue Dilator
- 1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW
- 1: Syringe: 5 mL Luer-Slip
- 1: Obturator: 8 Fr.
- 1: Cath-Gard® Catheter Contamination Shield: 80 cm with Sterile Tape
- 1: Dust Cap: Non-Vented
- 1: Stopcock: 4-way High-Flow

Not made with natural rubber latex.
Fluid path components are non-pyrogenic.

- fr** Jeu d'introduction de gaines par voie percutanée avec valve hémostatique/port latéral intégrés pour emploi avec les cathéters 7 - 8 Fr.
- de** Perkutanes Schleuseneinführbesteck mit eingebauter Hämostaseventil/ Seitenanschluß-Baugruppe zum Gebrauch mit Kathetern der Größe 7 - 8 Fr.
- it** Set per introduzione di guaina percutanea con valvola emostatica/bocchetta laterale incorporate per l'uso con cateteri da 7 - 8 Fr.
- pt** Conjunto de Introdução de Bainha Percutânea com Válvula de Hemostase/Integral/Porta Lateral, para utilizar com Cateteres 7 - 8 Fr.
- es** Conjunto introductor de vaina percutánea con válvula hemostática y orificio lateral integral para uso con catéteres de 7 - 8 Fr.
- sv** Set för perkutan skyddsintroducerare med integral hemostasventil/sidoöppning för bruk med katetrar mellan 7 - 8 Fr.



LOT LotNumber
YYYY-MM-DD



(01)XXXXXXXXXXXXXXXX
(17)YYMMDD
(10)LotNumber

Arrow International, Inc.
Subsidiary of Teleflex Incorporated
2400 Bernville Road
Reading, PA 19605 USA

EC REP Teleflex Medical
IDA Business and
Technology Park
Dublin Road, Athlone
0086 Co. Westmeath, Ireland

Teleflex®
www.teleflex.com
LBL011299 (06/18)



REF AD-09903-S

9 Fr. | **10** cm sheath length | **3.9** mm sheath O.D.

Percutaneous Sheath Introducer Set with Integral Hemostasis Valve/Side Port for use with 7 - 8 Fr. Catheters

