

Special Order Product



REF ES-04306-N (IPNXXXXXX)

16 Ga. | **1** Lumen | **16** cm catheter length | **.032** inch dia. spring-wire guide

CENTRAL VENOUS CATHETERIZATION SET with Blue FlexTip® Catheter

Contents:

- 1: Indwelling Catheter: 16 Ga. x 16 cm Radiopaque Polyurethane with Blue FlexTip®, Integral Extension Line, Extension Line Clamp, Integral Suture Wing
- 1: Spring-Wire Guide: .032" (0.81 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End – "J" Tip on Other)
- 1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW
- 1: Syringe: 5 mL Luer-Slip
- 1: Tissue Dilator: 7 Fr. (2.3 mm) x 10.2 cm
- 1: SecondSite™ Adjustable Hub: Fastener
- 1: SecondSite™ Adjustable Hub: Catheter Clamp


- fr** Jeu de cathétérisme veineux central avec cathéter Blue FlexTip
- de** Besteck für zentralen Venenkatheter mit Blue-FlexTip-Katheter
- it** Set per cateterismo della vena centrale con catetere Blue FlexTip
- pt** Conjunto de Cateterização Venosa Central com Cateter Blue FlexTip
- es** Conjunto de cateterización venosa central con catéter Blue FlexTip
- sv** Set för central venkatetrisering med Blue FlexTip-kateter

All components are CE 0086 unless otherwise noted.

Rx only. Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.

Fluid path components are non-pyrogenic.

Not made with natural rubber latex.



Lumen	Priming Volume* (ml)	Flow Rate† (ml/hr)
Distal (16 Ga.)	0.40	3700

* Priming volumes are approximate and are done without the injection cap. Injection cap priming volume is 0.17 ml.

† Flow rates are done with normal saline, room temperature, 100 cm head height and represent approximate flow capabilities.



LOT LotNumber
YYYY-MM-DD



(01)XXXXXXXXXXXXXX
(17)YYMMDD
(10)LotNumber

EC REP Teleflex Medical
IDA Business and
Technology Park
Dublin Road, Athlone
Co. Westmeath, Ireland

Arrow International, Inc.
Subsidiary of Teleflex Incorporated
2400 Bernville Road
Reading, PA 19605 USA

Teleflex
www.teleflex.com

LBL012057 (08/18)



REF ES-04306-N **16** Ga. | **1** Lumen | **16** cm catheter length | **.032** inch dia. spring-wire guide

CENTRAL VENOUS CATHETERIZATION SET with Blue FlexTip® Catheter

