

Special Order Product



REF CI-09800-SB (IPNXXXXXX)

8.5 Fr. | 10 cm sheath length | .035 inch dia. spring-wire guide

**PERCUTANEOUS SHEATH INTRODUCER SET with Integral Hemostasis Valve/Side Port for use with 7 - 7.5 Fr. Catheters**

Contents:

- 1: Sheath: 8.5 Fr. x 4-1/8" (10 cm) Arrow-Flex® Radiopaque Polyurethane with Integral Hemostasis Valve/Side Port and Tissue Dilator
- 1: Spring-Wire Guide: .035" (0.89 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End – "J" Tip on Other)
- 1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW
- 1: Obturator: 8 Fr.

Sheath O.D.: 3.6 mm

**Rx only. Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.**

This product contains dry natural rubber.

Fluid path components are non-pyrogenic.



LOT LotNumber  
YYYY-MM-DD



(01)XXXXXXXXXXXXXX  
(17)YYMMDD  
(10)LotNumber

EC REP Teleflex Medical  
IDA Business and  
Technology Park  
Dublin Road, Athlone  
0086 Co. Westmeath, Ireland

Arrow International, Inc.  
Subsidiary of Teleflex Incorporated  
2400 Bernville Road  
Reading, PA 19605 USA

Teleflex®  
www.teleflex.com  
LBL012208 (09/18)



REF CI-09800-SB

8.5 Fr. | 10 cm sheath length | .035 inch dia. spring-wire guide

PERCUTANEOUS SHEATH INTRODUCER SET with Integral Hemostasis Valve/Side Port for use with 7 - 7.5 Fr. Catheters

