

Dimensions: Width: (10.8 cm) Length: (11.4 cm)

PBZ-05000-109A

ARROW®

REF AN-04318-EG (IPNXXXXXX)

Introducer Needle: Echogenic 18 Ga. x 2-1/2" (6.35 cm) XTW

Contents: _____

1: Introducer Needle: Echogenic 18 Ga. x 2-1/2"
(6.35 cm) XTW

Do not use if package is damaged.

Not made with natural rubber latex.
Fluid path components are non-pyrogenic.

CE **REP** Teleflex Medical
IDA Business and
Technology Park
Dublin Road, Athlone
Co. Westmeath, Ireland
0086

Arrow International, Inc.
Subsidiary of Teleflex Incorporated
2400 Bernville Road
Reading, PA 19605 USA

Teleflex®
www.teleflex.com
LBL013067 (04/19)



LOT LotNumber
YYYY-MM-DD

(01)XXXXXXXXXXXXXXXXXX
(17)YYMMDD
(10)LotNumber



REF AN-04318-EG
LOT LotNumber YYYY-MM-DD
(01)XXXXXXXXXXXXXXXXXX

