



REF CS-45854-E (IPNXXXXXX)

8.5 Fr. 4 Lumen 20 cm catheter length .032 inch dia. spring-wire guide

Arrowg+ard Blue Plus® Four-Lumen CVC

Contents:

- 1: Arrowg+ard Blue Plus®¹ Four-Lumen Catheter: 8.5 Fr. (3.0 mm OD) x 20 cm
- 1: Spring-Wire Guide, Marked: .032" (0.81 mm) dia. x 23-5/8" (60 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer with ECG Mark
- 1: Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle
- 1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW
- 1: Pressure Transduction Probe
- 1: Arrow® Raulerson Spring-Wire Introduction Syringe: 5 mL
- 1: Tissue Dilator: 9 Fr. (3.0 mm) x 10.2 cm
- 4: Dust Cap: Non-Vented
- 1: SecondSite™ Adjustable Hub: Fastener
- 1: SecondSite™ Adjustable Hub: Catheter Clamp

- CS** Čtyřlumenový CVC Arrowg+ard Blue Plus
- nl** Arrowg+ard Blue Plus-CVC met vier lumina
- fr** CVC quadruple lumière Arrowg+ard Blue Plus
- de** Arrowg+ard Blue Plus vierlumiger ZVK
- it** CVC a quattro lumi Arrowg+ard Blue Plus
- pt** CVC de quatro lúmenes Arrowg+ard Blue Plus
- es** CVC de cuatro luces Arrowg+ard Blue Plus
- sv** Arrowg+ard Blue Plus fyrlumen-CVK

¹Licensed under US Patent No. 7,329,412.

All components are CE 0086 unless otherwise noted.

Rx only. Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.

California Prop. 65

⚠ **WARNING:** Cancer and Reproductive Harm.
www.P65Warnings.ca.gov

Not made with natural rubber latex. Store below 25°C (77°F). Avoid excessive heat above 40°C (104°F).

Fluid path components are non-pyrogenic.

Contraindications: The Arrowg+ard Blue Plus antimicrobial catheter is contraindicated for patients with known hypersensitivity to chlorhexidine, silver sulfadiazine and/or sulfa drugs.



Lumen	Priming Volume* (mL)	Flow Rate† (mL/hr)
Distal (16 Ga.)	0.40	2700
Medial 1 (14 Ga.)	0.50	5100
Medial 2 (18 Ga.)	0.36	1400
Proximal (18 Ga.)	0.39	1600

* Priming volumes are approximate and are done without accessories.

† Based on 100 cm head height.



LOT LotNumber

YYYY-MM-DD



(01)XXXXXXXXXXXXXX
(17)YYMMDD
(10)LotNumber



Teleflex Medical
IDA Business and
Technology Park
Dublin Road, Athlone
Co. Westmeath, Ireland

Arrow International, Inc.
Subsidiary of Teleflex Incorporated
2400 Bernville Road
Reading, PA 19605 USA

Teleflex®
www.teleflex.com

LBL015821 (08/20)



REF CS-45854-E 8.5 Fr. 4 Lumen 20 cm catheter length .032 inch dia. spring-wire guide

Arrowg+ard Blue Plus® Four-Lumen CVC

