



REF AK-09601

6 Fr. sheath

10 cm sheath length

.035 inch dia. spring-wire guide

Percutaneous Sheath Introducer Kit with Integral Hemostasis Valve/Side Port for use with 4 - 6 Fr. Catheters

Contents:

- 1: Sheath: 6 Fr. x 4" (10 cm) Arrow-Flex® Radiopaque Polyurethane with Integral Hemostasis Valve/Side Port and Tissue Dilator
- 1: Spring-Wire Guide: .035" (0.89 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other)
- 1: Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle
- 1: Injection Needle: 22 Ga. x 1-1/2" (3.81 cm)
- 1: Injection Needle: 25 Ga. x 1" (2.54 cm)
- 1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW
- 1: Syringe: 3 mL Luer-Lock
- 2: Syringe: 5 mL Luer-Slip
- 1: Obturator: 8 Fr.
- 1: Cath-Gard® Catheter Contamination Shield: 30.5 cm with Sterile Tape
- 1: SharpsAway® Disposal Cup
- 1: Drape: 24" x 36" (61 cm x 91 cm) with 4" (10 cm) fenestration
- 1: Scalpel: #11
- 2: Prep Sponge Swab
- 2: Gauze Pad: 2" x 2" (5 cm x 5 cm)
- 5: Gauze Pad: 4" x 4" (10 cm x 10 cm)
- 1: Suture: 3-0 Braided Silk with Straight Needle
- 1: Tray: Prep

Sheath O.D.: 2.8 mm

Rx only

Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.

Caution: This Product Contains Dry Natural Rubber.

Fluid path components are non-pyrogenic.

- fr** Kit d'introduction de gaines par voie percutanée avec valve hémostatique/port latéral intégrés pour emploi avec les cathéters 4 - 6 Fr.
- de** Perkutane Schleuseneinführungs- und Baugruppe zum Gebrauch mit Kathetern der Größen 4 - 6 Fr.
- it** Kit per introduzione di guaina percutanea con valvola emostatica/bocchetta laterale incorporate per l'uso con cateteri da 4 - 6 Fr.
- pt** Kit de Introdução de Bainha Percutânea com Dispositivo de Válvula de Hemostase Integral e Porta Lateral, para utilizar com Cateteres 4 - 6 Fr.
- es** Juego introductor de vaina percutánea con válvula hemostática y orificio lateral integrales para uso con catéteres de 4 a 6 Fr.
- sv** Uppsättning för perkutan skyddsintroducerare med integral hemostasventil/sidoöppning för bruk med katetrar på 4 och 6 Fr.



LOT



EC REP Teleflex Medical
IDA Business and
Technology Park
Dublin Road, Athlone
0086 Co. Westmeath, Ireland

Arrow International, Inc.
Subsidiary of Teleflex Incorporated
2400 Bernville Road
Reading, PA 19605 USA

Teleflex
www.teleflex.com

LBL015876 R01 (2021-05)



REF AK-09601

6 Fr. sheath

10 cm sheath length

.035 inch dia. spring-wire guide

Percutaneous Sheath Introducer Kit with Integral Hemostasis Valve/Side Port for use with 4 - 6 Fr. Catheters