



REF AK-09802-A

8.5 Fr. **10** cm sheath length **.035** inch dia. spring-wire guide

Percutaneous Sheath Introducer Kit with Integral Hemostasis Valve/Side Port for use with 7 - 7.5 Fr. Catheters

Contents:

- 1: Sheath: 8.5 Fr. x 4" (10 cm) Arrow-Flex® Radiopaque Polyurethane with Integral Hemostasis Valve/ Side Port and Tissue Dilator
- 1: Spring-Wire Guide: .035" (0.89 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other)
- 1: Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle
- 1: Injection Needle: 25 Ga. x 1" (2.54 cm) and 3 mL Luer-Lock Syringe
- 1: Injection Needle: 22 Ga. x 1-1/2" (3.81 cm) and 5 mL Luer-Slip Syringe
- 1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW and 5 mL Luer-Slip Syringe
- 1: Syringe: 1 mL Luer-Slip
- 1: 3 mL Applicator 2% CHG and 70% IPA Chloraprep®¹ One-Step Solution with Hi-Lite Orange™ Tint (non-sterile solution)
- 1: Obturator: 8 Fr.
- 1: Single-Lumen Infusion Catheter, SLIC®: 7 Fr. x 6" (16 cm) Radiopaque Polyurethane with Obturator
- 1: Cath-Gard® Catheter Contamination Shield: 80 cm with Sterile Tape
- 1: SharpsAway® Disposal Cup
- 3: Drape: 26" x 18" (66 cm x 46 cm) with adhesive
- 1: Scalpel: #11
- 1: Stopcock: 3-way High-Flow
- 5: Gauze Pad: 4" x 4" (10 cm x 10 cm)
- 1: Suture: 3-0 Nylon with Straight Needle

¹A registered trademark of CareFusion or one of its subsidiaries.

Rx only.

Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.

California Prop. 65

⚠ WARNING: Cancer and Reproductive Harm.
www.P65Warnings.ca.gov

This product contains dry natural rubber.

Store between 20 – 25°C (68 – 77°F). Avoid freezing and excessive heat above 40°C (104°F).

Sterilized using ethylene oxide unless otherwise indicated in the contents list.

Fluid path components are non-pyrogenic.



Arrow International, Inc.
Subsidiary of Teleflex Incorporated
2400 Bernville Road
Reading, PA 19605 USA

Teleflex®
www.teleflex.com
LBL015880 (05/19)



REF AK-09802-A

8.5 Fr. **10** cm sheath length **.035** inch dia. spring-wire guide