



**REF CV-13122-F**

(IPNXXXXXX)

**12** Fr. **2** Lumen **13** cm catheter length **.035** inch dia. spring-wire guide

## You-Bend™ Two-Lumen Hemodialysis Catheterization Set with Blue FlexTip® Catheter for High Volume Infusions

**Contents:** \_\_\_\_\_

- 1: Two-Lumen Indwelling Catheter‡: 12 Fr. x 13 cm Radiopaque Polyurethane with Blue FlexTip®, Extension Line Clamps
- 1: Spring-Wire Guide, Marked: .035" (0.89 mm) dia. x 23-5/8" (60 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer
- 1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW
- 1: Syringe: 5 mL Luer-Slip
- 1: Tissue Dilator: 12 Fr. (4.0 mm) x 14 cm
- 2: Dust Cap: Non-Vented


**Rx only.**

**Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.**

**‡MR Conditional**

Not made with natural rubber latex.

Fluid path components are non-pyrogenic.



Lumen	Priming Volume* (mL)	Flow Rate† (mL/min)	Venous Pressure Setting (mmHg)
Distal (12 Ga.)	1.5	230	100
		400	200
Proximal (12 Ga.)	1.4	230	100
		400	200

\* Priming volumes are approximate and are done without accessories.

† Blood flow rates are approximate and were determined by varying the venous pressure setting from 100 to 200 mmHg in an in vitro experiment using blood analog solution.



Arrow International, Inc.  
Subsidiary of Teleflex Incorporated  
2400 Bernville Road  
Reading, PA 19605 USA

**Teleflex®**  
www.teleflex.com

LBL015960 (07/19)

**LOT** LotNumber

YYYY-MM-DD



(01)XXXXXXXXXXXXXX  
(17)YYMMDD  
(10)LotNumber



**REF CV-13122-F 12** Fr. **2** Lumen **13** cm catheter length **.035** inch dia. spring-wire guide

