

Dimensions:
Width: 6-9/64" (15.6 cm)
Length: 8-21/32" (22 cm)

PB-50014-100B
PBZ-50014-103A

T50014100A6 Dimensions: 5 7/8" x 7 5/8"

ARROW®



REF CV-15854 (IPNXXXXXX)

8.5 Fr. 4 Lumen 20 cm catheter length .032 inch dia. spring-wire guide

Quad-Lumen Central Venous Catheterization Set with Blue FlexTip® Catheter

Contents:

- 1: Quad-Lumen Indwelling Catheter: 8.5 Fr. x 20 cm Radiopaque Polyurethane with Blue FlexTip®, Extension Line Clamps
- 1: Spring-Wire Guide, Marked: .032" (0.81 mm) dia. x 23-5/8" (60 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer with ECG Mark
- 1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW
- 1: Syringe: 5 mL Luer-Slip
- 1: Tissue Dilator: 9 Fr. (3.0 mm) x 10.2 cm
- 4: Dust Cap: Non-Vented
- 1: SecondSite™ Adjustable Hub: Fastener
- 1: SecondSite™ Adjustable Hub: Catheter Clamp

All components are CE 0086 unless otherwise noted.

Rx only. **Warning:** Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.

California Prop. 65


⚠ **WARNING:** Cancer and Reproductive Harm.

www.P65Warnings.ca.gov

Not made with natural rubber latex.

Fluid path components are non-pyrogenic.

- fr** Jeu de cathétérisme veineux central à quatre lumières avec cathéter Blue FlexTip
- de** Besteck für vierlumigen, zentralen Venenkatheter mit Blue FlexTip Katheter
- it** Set per cateterismo della vena centrale a quattro lumi con catetere Blue FlexTip
- pl** Zestaw do cewnikowania żyły centralnej z cewnikiem 4-kanalowym z końcówką Blue FlexTip
- pt** Conjunto de Cateterização Venosa Central com Lúmen Quádruplo e Cateter Blue FlexTip
- ru** Четырех-просветный комплект для установки центрального венозного катетера с катетером Blue FlexTip
- sl** Štiri-lumenski centralni venski komplet za kateterizacijo s katetrom Blue FlexTip
- es** Conjunto de cateterización venosa central de cuatro luces con catéter Blue FlexTip
- sv** Set för fyrkanalig central venkatetrering med Blue FlexTip kateter
- tr** Dört Lümen Santral Venöz Kateterizasyon Seti, Mavi esnek uçlu Kateterli



Catheter Cross Section

Lumen	Priming Volume* (mL)	Flow Rate† (mL/hr)
Distal (16 Ga.)	0.43	2700
Medial 1 (14 Ga.)	0.54	5100
Medial 2 (18 Ga.)	0.36	1400
Proximal (18 Ga.)	0.41	1600

* Priming volumes are approximate and are done without accessories.

† Flow rates are done with normal saline, room temperature, 100 cm head height and represent approximate flow capabilities.



LOT LotNumber
YYYY-MM-DD



(01)XXXXXXXXXXXXXX
(17)YYMMDD
(10)LotNumber

EC REP Teleflex Medical
IDA Business and
Technology Park
Dublin Road, Athlone
Co. Westmeath, Ireland
0086

Arrow International, Inc.
Subsidiary of Teleflex Incorporated
2400 Bernville Road
Reading, PA 19605 USA

teleflex®
www.teleflex.com

LBL015963 (07/19)



8.5 Fr. 4 Lumen 20 cm catheter length .032 inch dia. spring-wire guide

REF CV-15854

Quad-Lumen Central Venous Catheterization Set with Blue FlexTip® Catheter

ARROW®