

Dimensions:
Width: 6-9/64" (15.6 cm)
Length: 8-21/32" (22 cm)

PB-50014-100B
PBZ-50014-103A

T50014100A6 Dimensions: 5 7/8" x 7 5/8"

ARROW®



REF CV-25854 (IPNXXXXXX)

8.5 Fr. 4 Lumen 20 cm catheter length .032 inch dia. spring-wire guide

Quad-Lumen Central Venous Catheterization Set with Blue FlexTip® ARROWg+ard Blue® Catheter

Contents:

- 1: Quad-Lumen Indwelling Catheter: 8.5 Fr. x 20 cm Radiopaque Polyurethane with Blue FlexTip®, ARROWg+ard® Antimicrobial Surface Treatment¹, Extension Line Clamps
- 1: Spring-Wire Guide, Marked: .032" (0.81 mm) dia. x 23-5/8" (60 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer with ECG Mark
- 1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW
- 1: Syringe: 5 mL Luer-Slip
- 1: Tissue Dilator: 9 Fr. (3.0 mm) x 10.2 cm
- 4: Dust Cap: Non-Vented
- 1: SecondSite™ Adjustable Hub: Fastener
- 1: SecondSite™ Adjustable Hub: Catheter Clamp

¹Licensed under US Patent No. 6,706,024.

All components are CE 0086 unless otherwise noted.

- fr** Jeu de cathétérisme veineux central à quatre lumières avec cathéter Blue FlexTip ARROWg+ard Blue
- de** Besteck für vierlumigen, zentralen Venenkatheter mit Blue FlexTip ARROWg+ard Blue Katheter
- it** Set per cateterismo della vena centrale a quattro lumi con catetere Blue FlexTip ARROWg+ard Blue
- pl** Zestaw do cewnikowania żyły centralnej z 4-kanalowym cewnikiem z końcówką Blue FlexTip i powłoką antybakteryjną ARROWg+ard Blue
- pt** Conjunto de Cateterização Venosa Central com Lúmen Quádruplo e Cateter Blue FlexTip, ARROWg+ard Blue
- ru** Четырех-просветный набор для центральной венозной катетеризации с катетером Blue FlexTip ARROWg+ard Blue
- sl** Štiri-lumenski centralni venski komplet za katetrizacijo s katetrom Blue FlexTip ARROWg+ard Blue
- es** Conjunto de cateterización venosa central de cuatro luces con catéter Blue FlexTip ARROWg+ard Blue
- sv** Set för fyrkanalig central venkatetrisering med Blue FlexTip ARROWg+ard Blue kateter
- tr** Dört Lümenli Santral Venöz Kateterizasyon Seti, Mavi esnek uçlu ARROWg+ard Blue Kateterli

Rx only. **Warning:** Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.


California Prop. 65

⚠ **WARNING:** Cancer and Reproductive Harm.
www.P65Warnings.ca.gov

Not made with natural rubber latex. Store below 25°C (77°F). Avoid excessive heat above 40°C (104°F).

Fluid path components are non-pyrogenic.

Contraindications: The Arrowg+ard Blue antimicrobial catheter is contraindicated for patients with known hypersensitivity to chlorhexidine, silver sulfadiazine and/or sulfa drugs (refer to product instructions for references).



Lumen	Priming Volume* (mL)	Flow Rate† (mL/hr)
Distal (16 Ga.)	0.43	2700
Medial 1 (14 Ga.)	0.54	5100
Medial 2 (18 Ga.)	0.36	1400
Proximal (18 Ga.)	0.41	1600

* Priming volumes are approximate and are done without accessories.

† Based on 100 cm head height.



LOT LotNumber
YYYY-MM-DD



(01)XXXXXXXXXXXXXX
(17)YYMMDD
(10)LotNumber

EC REP Teleflex Medical
IDA Business and
Technology Park
Dublin Road, Athlone
Co. Westmeath, Ireland
0086

Arrow International, Inc.
Subsidiary of Teleflex Incorporated
2400 Bernville Road
Reading, PA 19605 USA

teleflex®
www.teleflex.com

LBL015966 (07/19)



REF CV-25854

8.5 Fr. 4 Lumen 20 cm catheter length .032 inch dia. spring-wire guide

Quad-Lumen Central Venous Catheterization Set with Blue FlexTip® ARROWg+ard Blue® Catheter

ARROW®