



REF SI-09875-E

ARROW®

8.5 Fr. 10 cm sheath length .035 inch dia. spring-wire guide

## Percutaneous Sheath Introducer Set with Integral Hemostasis Valve/ Side Port for use with 7.5 Fr. Catheters

### Contents:

- 1: Sheath: 8.5 Fr. x 4" (10 cm) Arrow-Flex® Radiopaque Polyurethane with Integral Hemostasis Valve/Side Port and Tissue Dilator
- 1: Spring-Wire Guide, Marked: .035" (0.89 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer
- 1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW
- 1: Pressure Transduction Probe
- 1: Arrow® Raulerson Spring-Wire Introduction Syringe: 5 mL
- 1: Obturator: 8 Fr.
- 1: Cath-Gard® Catheter Contamination Shield: 80 cm with TwistLock™ Adapter
- 3: Gauze Pad: 4" x 4" (10 cm x 10 cm)

Sheath O.D.: 3.6 mm

Rx only.

**Warning:** Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.

This Product Contains Dry Natural Rubber.

Fluid path components are non-pyrogenic.

- cs** Set s perkutánním zaváděcím pouzdem a integrovaným hemostatickým ventilem a bočním portem pro použití s katetry 7,5 Fr.
- da** Perkutan sheath introducer-sæt med integreret hæmostaseventil/ sideåbning til brug med 7,5 Fr. katetre
- nl** Percutane hulsintroducerset met integrale hemostaseklep/zijpoort voor gebruik met 7,5 Fr. katheters
- fr** Jeu d'introduction de gaines par voie percutanée avec valve hémostatique/ port latéral intégrés pour emploi avec les cathéters 7,5 Fr.
- de** Perkutanes Schleuseneinführbesteck mit eingebauter Hämostaseventil/ Seitenanschluß- Baugruppe zum Gebrauch mit Kathetern der Größen 7,5 Fr.
- el** Σετ διαδερμικού εισαγωγέα θηκαριού με ενιαία αιμοστατική βαλβίδα/ πλευρική θύρα για χρήση με καθετήρες 7,5 Fr.
- hu** Percutan hüvelybevezető felszerelés egybeépített vérzéscsillapító szeleppel és oldalcsatlakozóval 7,5 Fr. katéterekkel való használatra
- it** Set per introduzione di guaina percutanea con valvola emostatica/ bocchetta laterale incorporate per l'uso con cateteri da 7,5 Fr.
- pl** Zestaw do przezskórnego wprowadzania koszulki, z integralnym zaworem hemostaticznym/portem bocznym do użycia z cewnikami o rozmiarze 7,5 Fr.
- pt** Conjunto de Introdução de Bainha Percutânea com Dispositivo de Válvula de Hemostase Integral e Porta Lateral, para utilizar com Cateteres 7,5 Fr.
- ru** Комплект чрескожной капсулы-интубатора с несъемным гемостазным клапаном/боковым портом, для использования с катетерами 7,5 Fr.
- sk** Perkutánná sada zavádzača puzdra s integrálnym hemostatickým ventilom/ postranným portom pre použitie s 7,5 Fr. Katétrami
- sl** Uvajalni komplet s perkutanim tulcem z vgrajenim hemostatskim ventilom/ stranskim vhomom za uporabo s katetri 7,5 Fr.
- es** Conjunto introductor de vaina percutánea con válvula hemostática y orificio lateral integrales para uso con catéteres de 7,5 Fr.
- sv** Set för perkutan skyddsintroducerare med integral hemostasventil/ sidoöppning för bruk med katetrar mellan 7,5 Fr.
- tr** Perkütan Kılıf İntroduser Seti, Entegre Hemostaz Valfi/Yan Portlu, 7,5 Fr. Kateterler ile kullanılmak üzere



Arrow International, Inc.  
Subsidiary of Teleflex Incorporated  
2400 Bernville Road  
Reading, PA 19605 USA



**Teleflex®**  
www.teleflex.com

LBL016048 (09/19)



REF SI-09875-E

8.5 Fr. 10 cm sheath length .035 inch dia. spring-wire guide

## Percutaneous Sheath Introducer Set with Integral Hemostasis Valve/Side Port for use with 7.5 Fr. Catheters

ARROW®