



REF SI-09880

ARROW®

8.5 Fr. 10 cm sheath length .035 inch dia. spring-wire guide

Percutaneous Sheath Introducer Set with Integral Hemostasis Valve/ Side Port for use with 7 - 7.5 Fr. Catheters

Contents:

- 1: Sheath: 8.5 Fr. x 4" (10 cm) Arrow-Flex® Radiopaque Polyurethane with Integral Hemostasis Valve/Side Port and Tissue Dilator
- 1: Spring-Wire Guide: .035" (0.89 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other)
- 1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW
- 1: Obturator: 8 Fr.
- 1: Cath-Gard® Catheter Contamination Shield: 80 cm with Sterile Tape
- 3: Gauze Pad: 4" x 4" (10 cm x 10 cm)

Rx only.

Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.

Not made with natural rubber latex.

Fluid path components are non-pyrogenic.

- fr** Jeu d'introduction de gaines par voie percutanée avec valve hémostatique/port latéral intégrés pour emploi avec les cathéters 7 - 7,5 Fr.
- de** Perkutanes Schleuseneinführbesteck mit eingebauter Hämostaseventil/Seitenanschluß-Baugruppe zum Gebrauch mit Kathetern der Größen 7 - 7,5 Fr.
- it** Set per introduzione di guaina percutanea con valvola emostatica/ bocchetta laterale incorporate per l'uso con cateteri da 7 - 7,5 Fr.
- pl** Zestaw do przezskórnego wprowadzania koszulki, z integralnym zaworem hemostatycznym / portem bocznym, do użycia z cewnikami o rozmiarze 7 - 7,5 Fr.
- pt** Conjunto de Introdução de Bainha Percutânea com Dispositivo de Válvula de Hemostase Integral e Porta Lateral, para utilizar com Cateteres 7 - 7,5 Fr.
- ru** Комплект чрескожной капсулы-интубатора с несъемным гемостазным клапаном/боковым портом, для использования с катетерами 7 - 7,5 Fr.
- sl** Uvajalni komplet s perkutanim tulcem z vgrajenim hemostatskim ventilom/stranskim vhodom za uporabo s katetri 7 – 7,5 Fr.
- es** Conjunto introductor de vaina percutánea con válvula hemostática y orificio lateral integrales para uso con catéteres de 7 a 7,5 Fr.
- sv** Set för perkutan skyddsintroducerare med integral hemostasventil/ sidoöppning för bruk med katetrar mellan 7 och 7,5 Fr.
- tr** Perkütan Kılıf Introducer Seti, Entegre Hemostaz Valfi/Yan Portlu, 7 - 7,5 Fr. Kateterler ile kullanılmak üzere



Arrow International, Inc.
Subsidiary of Teleflex Incorporated
2400 Bernville Road
Reading, PA 19605 USA



Teleflex®
www.teleflex.com

LBL016050 (09/19)



REF SI-09880

8.5 Fr. 10 cm sheath length .035 inch dia. spring-wire guide

Percutaneous Sheath Introducer Set with Integral Hemostasis Valve/Side Port for use with 7 - 7.5 Fr. Catheters

ARROW®