



REF **CV-15123-F**

12 Fr. **3** Lumen **20** cm catheter length **.035** inch dia. spring-wire guide

Arrow-Howes™ Large-Bore Multi-Lumen Central Venous Catheterization Set with Blue FlexTip® Catheter for High Volume Infusions

Contents:

- 1: Multi-Lumen Indwelling Catheter: 12 Fr. x 20 cm Radiopaque Polyurethane with Blue FlexTip®, Extension Line Clamps
- 1: Spring-Wire Guide, Marked: .035" (0.89 mm) dia. x 26-3/4" (68 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer
- 1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW and 5 mL Luer-Slip Syringe
- 1: Tissue Dilator: 12 Fr. (4.0 mm) x 14 cm
- 3: Dust Cap: Non-Vented

All components are CE 0086 unless otherwise noted.


Rx only.

Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.

Not made with natural rubber latex.

Fluid path components are non-pyrogenic.

- fr** Jeu de cathétérisme veineux central de grand calibre à lumières multiples Arrow-Howes avec cathéter Blue FlexTip pour les perfusions de grand volume
- de** Arrow-Howes-Besteck für viellumigen, zentralen, großkalibrigen Venenkatheter mit Blue-FlexTip-Katheter für Infusionen großer Flüssigkeitsmengen
- it** Set per cateterismo della vena centrale multilume di diametro ampio Arrow-Howes con catetere Blue FlexTip per volumi di infusione elevati
- pt** Conjunto de Cateterização Venosa Central Arrow-Howes com Lúmen Múltiplo e Grande Calibre, com Cateter Blue FlexTip, para Infusões de Grandes Volumes
- es** Conjunto de cateterización venosa central de luz múltiple y gran diámetro Arrow-Howes con catéter Blue FlexTip para infusiones de gran volumen
- sv** Arrow-Howes set för storhålig flerkanalig central venkatetrering med Blue FlexTip-kateter för infusioner med hög volym



Lumen	Priming Volume* (mL)	Flow Rate† (mL/hr)
Distal (16 Ga.)	0.4	80
Medial (12 Ga.)	1.6	310
Proximal (12 Ga.)	1.5	330

* Priming volumes are approximate and are done without accessories.
† Based on 100 cm head height.



Arrow International, Inc.
Subsidiary of Teleflex Incorporated
2400 Bernville Road
Reading, PA 19605 USA



REF **CV-15123-F**

Arrow-Howes™ Large-Bore Multi-Lumen Central Venous Catheterization Set with Blue FlexTip® Catheter for High Volume Infusions

12 Fr. **3** Lumen **20** cm catheter length **.035** inch dia. spring-wire guide

