



REF CS-24301-E

16^{Ga.} 1^{Lumen} 20^{cm} catheter length .032^{inch dia. spring-wire guide}

Central Venous Catheterization Set with Blue FlexTip® ARROWg+ard Blue® Catheter

Contents:

- 1: Indwelling Catheter: 16 Ga. x 20 cm Radiopaque Polyurethane with Blue FlexTip®, ARROWg+ard® Antimicrobial Surface Treatment¹, Integral Extension Line, Extension Line Clamp, Integral Suture Wing
- 1: Spring-Wire Guide, Marked: .032" (0.81 mm) dia. x 23-5/8" (60 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer
- 1: Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle
- 1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW
- 1: Pressure Transduction Probe
- 1: Arrow® Raulerson Spring-Wire Introduction Syringe: 5 mL
- 1: Tissue Dilator: 7 Fr. (2.3 mm) x 10.2 cm
- 1: Dust Cap: Non-Vented
- 1: SecondSite™ Adjustable Hub: Fastener
- 1: SecondSite™ Adjustable Hub: Catheter Clamp

¹Licensed under US Patent No. 6,706,024.

All components are CE 0086 unless otherwise noted.


Rx only. Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.

Not made with natural rubber latex. Store below 25°C (77°F). Avoid excessive heat above 40°C (104°F).

Fluid path components are non-pyrogenic.

Contraindications: The Arrowg+ard Blue antimicrobial catheter is contraindicated for patients with known hypersensitivity to chlorhexidine, silver sulfadiazine and/or sulfa drugs (refer to product instructions for references).

- fr** Jeu de cathétérisme veineux central avec cathéter Blue FlexTip ARROWg+ard Blue
- de** Besteck für zentralen Venenkatheter mit Blue-FlexTip-ARROWg+ard-Blue-Katheter
- it** Set per cateterismo della vena centrale con catetere Blue FlexTip ARROWg+ard Blue
- pl** Zestaw do cewnikowania żył centralnych z cewnikiem Blue FlexTip ARROWg+ard Blue
- pt** Conjunto de Cateterização Venosa Central com Cateter Blue FlexTip, ARROWg+ard Blue
- ru** Комплект для катетеризации главных вен с катетером Blue FlexTip ARROWg+ard Blue
- sl** Komplet za centralno vensko katetrizacijo s katetrom Blue FlexTip ARROWg+ard Blue
- es** Conjunto para cateterización venosa central con catéter Blue FlexTip ARROWg+ard Blue
- sv** Set för central venkatetrisering med Blue FlexTip ARROWg+ard Blue-kateter
- tr** Santral Venöz Kateterizasyon Seti, Blue FlexTip ARROWg+ard Blue Kateter ile



Lumen	Priming Volume* (mL)	Flow Rate† (mL/hr)
Distal (16 Ga.)	0.43	3650

* Priming volumes are approximate and are done without accessories.

† Flow rates are done with normal saline, room temperature, 100 cm head height and represent approximate flow capabilities.



Arrow International, Inc.
Subsidiary of Teleflex Incorporated
2400 Bernville Road
Reading, PA 19605 USA

Teleflex®
www.teleflex.com

LBL026711 (06/19)



REF CS-24301-E

16^{Ga.} 1^{Lumen} 20^{cm} catheter length .032^{inch dia. spring-wire guide}

Central Venous Catheterization Set with Blue
FlexTip® ARROWg+ard Blue® Catheter

