



**REF CV-15122-FBZ**

**12** Fr. **2** Lumen **20** cm catheter length **.035** inch dia. spring-wire guide

## Two-Lumen Hemodialysis Catheterization Set with Blue FlexTip® Catheter for High Volume Infusions

**Contents:**

- 1: Two-Lumen Indwelling Catheter: 12 Fr. x 20 cm Radiopaque Polyurethane with Blue FlexTip®, Extension Line Clamps
- 1: Spring-Wire Guide, Marked: .035" (0.89 mm) dia. x 26-3/4" (68 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer
- 1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW
- 1: Syringe: 5 mL Luer-Slip
- 1: Tissue Dilator: 12 Fr. (4.0 mm) x 14 cm
- 2: Injection Site Cap

- pt** Conjunto de Cateterização para Hemodiálise com Lúmen Duplo e Cateter Blue FlexTip, para Infusões de Grandes Volumes
- es** Conjunto de cateterización de dos lúmenes para hemodiálisis con catéter Blue FlexTip para infusiones de gran volumen

All components are CE 0086 unless otherwise noted.

**Rx only.**

**Warning:** Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.


**California Prop. 65**

**⚠ WARNING:** Cancer and Reproductive Harm. [www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov)

Not made with natural rubber latex.

Fluid path components are non-pyrogenic.

**Contraindications:** The Arrow large-bore two-lumen catheter is not designed for long-term hemodialysis or for use in patients with thrombosed vessels.



Arrow International, Inc.  
 Subsidiary of Teleflex Incorporated  
 2400 Bernville Road  
 Reading, PA 19605 USA

**Teleflex®**  
[www.teleflex.com](http://www.teleflex.com)

LBL027217 (07/19)



**REF CV-15122-FBZ** **12** Fr. **2** Lumen **20** cm catheter length **.035** inch dia. spring-wire guide