



REF CV-45703 (IPNXXXXXX)

7 Fr. 3 Lumen 20 cm catheter length .032 inch dia. spring-wire guide

## ARROWg+ard Blue Plus® Multi-Lumen CVC Set

**Contents:**

- 1: Multi-Lumen Indwelling Catheter: 7 Fr. x 20 cm Radiopaque Polyurethane with Blue FlexTip®, ARROWg+ard Blue Plus® Antimicrobial Surface Treatment<sup>1</sup>, Extension Line Clamps
- 1: Spring-Wire Guide, Marked: .032" (0.81 mm) dia. x 23-5/8" (60 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer with ECG Mark
- 1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW
- 1: Syringe: 5 ml Luer-Slip
- 1: Tissue Dilator: 8.5 Fr. (2.8 mm) x 10.2 cm
- 3: Dust Cap: Non-Vented
- 1: SecondSite™ Adjustable Hub: Fastener
- 1: SecondSite™ Adjustable Hub: Catheter Clamp

<sup>1</sup>Licensed under US Patent Nos. 6,706,024, 6,872,195 and 7,329,412.

All components are CE 0086 unless otherwise noted.


- cs** Vícelumenový CVC set ARROWg+ard Blue PLUS
- nl** ARROWg+ard Blue PLUS multilumen centraal-veneuze katheterset
- fr** Jeu de cathétères veineux centraux à lumières multiples ARROWg+ard Blue PLUS
- de** Set mit mehrlumigen ARROWg+ard Blue PLUS ZVK
- it** Set CVC multilume ARROWg+ard Blue PLUS
- pt** Conjunto de CVC de Lúmen Múltiplo ARROWg+ard Blue PLUS
- es** Conjunto de CVC de luz múltiple ARROWg+ard Blue PLUS
- sv** ARROWg+ard Blue PLUS flerlumen CVK-set

**Rx only. Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.**

Not made with natural rubber latex. Store below 25°C (77°F). Avoid excessive heat above 40°C (104°F).

Fluid path components are non-pyrogenic.

**Contraindications: The Arrowg+ard Blue PLUS antimicrobial catheter is contraindicated for patients with known hypersensitivity to chlorhexidine, silver sulfadiazine and/or sulfa drugs (refer to product instructions for references).**



Lumen	Priming Volume* (mL)	Flow Rate† (mL/hr)
Distal (16 Ga.)	0.44	2283
Medial (18 Ga.)	0.39	913
Proximal (18 Ga.)	0.39	996

\* Priming volumes are approximate and are done without accessories.

† Flow rate values are approximate and are determined using deionized water at 100 cm head height.



**LOT** LotNumber  
YYYY-MM-DD



(01)XXXXXXXXXXXXXX  
(17)YYMMDD  
(10)LotNumber

**EC REP** Teleflex Medical  
IDA Business and  
Technology Park  
Dublin Road, Athlone  
Co. Westmeath, Ireland



Arrow International, Inc.  
Subsidiary of Teleflex Incorporated  
2400 Bernville Road  
Reading, PA 19605 USA

**Teleflex®**  
www.teleflex.com

LBL027338 (12/19)



REF CV-45703

7 Fr. 3 Lumen 20 cm catheter length .032 inch dia. spring-wire guide

ARROWg+ard Blue Plus® Multi-Lumen CVC Set

