



REF CS-18763-E

7 Fr. 3 Lumen 60 cm catheter length .035 inch dia. spring-wire guide

# Multi-Lumen Antecubital Central Venous Catheterization Set with Blue FlexTip® Catheter

**Contents:**

- 1: Multi-Lumen Indwelling Catheter: 7 Fr. x 60 cm Radiopaque Polyurethane with Blue FlexTip®, Extension Line Clamps
- 1: Spring-Wire Guide: .035" (0.89 mm) dia. x 39-1/2" (100 cm) PTFE-Coated (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer
- 1: Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle
- 1: Pressure Transduction Probe
- 1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW and 5 mL Arrow® Raulerson Spring-Wire Introduction Syringe
- 1: Tissue Dilator: 8.5 Fr. (2.8 mm) x 10.2 cm
- 1: Dust Cap: Non-Vented
- 3: Dust Cap: Non-Vented
- 1: SecondSite™ Adjustable Hub: Fastener
- 1: SecondSite™ Adjustable Hub: Catheter Clamp
- 1: Scalpel: #11
- 1: Stopcock: 3-way

All components are CE 0086 unless otherwise noted.

**Rx only. Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.**


**California Prop. 65**

⚠ **WARNING:** Cancer and Reproductive Harm. [www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov)

Not made with natural rubber latex.

Fluid path components are non-pyrogenic.

- fr** Jeu de cathétérisme veineux central antébrachial à lumières multiples avec cathéter Blue FlexTip
- de** Besteck für viellumigen, zentralen Venenkatheter zur Einführung in der Ellenbeuge mit Blue-FlexTip-Katheter
- it** Set per cateterismo della vena centrale cubitale anteriore multilume con catetere Blue FlexTip
- pl** Wielokanałowy zestaw do cewnikowania przedłokciowych żył centralnych z cewnikiem Blue FlexTip
- pt** Conjunto de Cateterização Venosa Central Antecubital com Lúmen Múltiplo e Cateter Blue FlexTip
- ru** Комплект для катетеризации главных вен с многопросветным катетером Blue FlexTip для введения спереди от локтевого сустава
- sl** Komplet za antekubitalno centralno venosko katetrizacijo z več svetlinami s katetrom Blue FlexTip
- es** Conjunto de cateterización venosa central antecubital de luz múltiple con catéter Blue FlexTip
- sv** Set för flerkanalig antekubital central venkatetrering med Blue FlexTip-kateter
- tr** Multi-Lümen Antekübital Santral Venöz Kateterizasyon Seti, Blue FlexTip Kateter ile



Lumen	Priming Volume* (mL)	Flow Rate† (mL/hr)
Distal (16 Ga.)	0.70	1300
Medial (18 Ga.)	0.57	500
Proximal (18 Ga.)	0.57	500

\* Priming volumes are approximate and are done without accessories.

† Flow rates are done with normal saline, room temperature, 100 cm head height and represent approximate flow capabilities.



LOT



CE REP Teleflex Medical  
IDA Business and  
Technology Park  
Dublin Road, Athlone  
Co. Westmeath, Ireland  
0086

Arrow International, Inc.  
Subsidiary of Teleflex Incorporated  
2400 Bernville Road  
Reading, PA 19605 USA

**Teleflex®**  
[www.teleflex.com](http://www.teleflex.com)  
LBL027971 (02/20)



REF CS-18763-E

7 Fr. 3 Lumen 60 cm catheter length .035 inch dia. spring-wire guide

# Multi-Lumen Antecubital Central Venous Catheterization Set with Blue FlexTip® Catheter