



REF ES-04701

14 Ga. **1** Lumen **20** cm catheter length **.032** inch dia. spring-wire guide

Central Venous Catheterization Set with Blue FlexTip® Catheter

Contents:

- 1: Indwelling Catheter: 14 Ga. x 20 cm Radiopaque Polyurethane with Blue FlexTip®, Integral Extension Line, Extension Line Clamp, Side Holes, Integral Suture Wing
- 1: Spring-Wire Guide, Marked: .032" (0.81 mm) dia. x 23-5/8" (60 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer
- 1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW
- 1: Pressure Transduction Probe
- 1: Arrow® Raulerson Spring-Wire Introduction Syringe: 5 mL
- 1: Tissue Dilator: 7 Fr. (2.3 mm) x 10.2 cm
- 1: SecondSite™ Adjustable Hub: Fastener
- 1: SecondSite™ Adjustable Hub: Catheter Clamp


- fr** Jeu de cathétérisme veineux central avec cathéter Blue FlexTip
- de** Besteck für zentralen Venenkatheter mit Blue-FlexTip-Katheter
- it** Set per cateterismo della vena centrale con catetere Blue FlexTip
- pl** Zestaw do cewnikowania żył centralnych z cewnikiem Blue FlexTip
- pt** Conjunto de Cateterização Venosa Central com Cateter Blue FlexTip
- ru** Комплект для катетеризации главных вен с катетером Blue FlexTip
- sl** Komplet za centralno vensko katetrizacijo s katetrom Blue FlexTip
- es** Conjunto de cateterización venosa central con catéter Blue FlexTip
- sv** Set för central venkatetrisering med Blue FlexTip-kateter
- tr** Santral Venöz Kateterizasyon Seti, Blue FlexTip Kateter ile

All components are CE 0086 unless otherwise noted.

Rx only. Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.

Not made with natural rubber latex.

Fluid path components are non-pyrogenic.



Lumen	Priming Volume* (mL)	Flow Rate† (mL/hr)
Distal (14 Ga.)	0.55	4900

* Priming volumes are approximate and are done without the injection cap. Injection cap priming volume is 0.17 mL.

† Flow rates are done with normal saline, room temperature, 100 cm head height and represent approximate flow capabilities.



LOT



EC REP Teleflex Medical
IDA Business and
Technology Park
Dublin Road, Athlone
Co. Westmeath, Ireland



Arrow International, Inc.
Subsidiary of Teleflex Incorporated
2400 Bernville Road
Reading, PA 19605 USA

Teleflex®
www.teleflex.com

LBL027975 (01/20)



REF ES-04701

14 Ga. **1** Lumen **20** cm catheter length **.032** inch dia. spring-wire guide

Central Venous Catheterization Set with Blue FlexTip® Catheter

